



التحاليل البكتيريولوجية

كيميائى مياه- الدرجة الثالثة

تم إعداد المادة بواسطة الشركة القابضة لمياه الشرب والصرف الصحي قطاع تنمية الموارد البشرية - الادارة العامة لتخطيط المسار الوظيفي الإصدار الثاني - 2019.

, is sere

ار الوظي

Г Г • 1× 10

Table of contents

	1. INTRODUCTION	5
	1.1. Definitions	5
	1.2. Why Water microbiology is important?	5
	1.3. History of Water Microbiology	6
	1.4. Bacteriological Indicators	6
	2. Methods of Analysis	9
	2.1. Estimation of Heterotrophic Bacteria by Pour Plate Technique	9
	2.2. Estimation of Heterotrophic Bacteria by Spread Plate Method	17
	2.3. Estimation of Coliform Group by Multiple Tube Fermentation Technique (For Drinking Water)	24
	2.4. Estimation of Coliform Group by Multiple Tube Fermentation (For water other than Drinking Water)	34
	2.5. Detection and Enumeration of Total Coliform by Membrane Filter Techni	<u>que</u> 47
	2.6 .Detection and Enumeration of Thermotolerant (Fecal) Coliform by Membrane Filter Technique	57
	2.7. Detection and Enumeration of Fecal Enterococcus /Streptococcus Groups Membrane Filter Technique	<u>by</u> 66
An	inexes	75

مقدمة الإصدار الثانى

تهدف مجموعة البرامج التدريبية المعدة من إدارة المسار الوظيفي بالشركة القابضة لمياه الشرب والصرف الصحى إلى رفع كفاءة الكيميائين العاملين بالشركة القابضة لمياه الشرب والصرف الصحى والشركات التابعة لها وتنمية مهاراتهم ومعارفهم بالشكل الذي يضمن الوصول إلى كوب مياه نظيف وبيئة آمنة يرضى متطلبات وإحتياجات العملاء الكرام.

ويعتبر الإصدار الثانى من برامج المسار الوظيفى لوظيفة كيميائى مياه الشرب هو ثمرة جهود الكيميائيين العاملين بمعامل الشركات التابعة والمعمل المرجعى لمياه الشرب بالشركة القابضة بما تحمله من مزيج متجانس من الخبرات والكفاءات الذين لم يدخروا جهدا حتى يخرج هذا العمل بالطريقة اللائقة.

وجدير بالذكر أن هذا الإصـدار يعتبر مكتبة مرجعية وافية وشـاملة لجميع الجدارات المتضـمنة المهارات والمعارف التي تجعل الكيميائي كفؤا لوظيفته.

ومما تجدر الإشارة إليه بأنه تم الاعتماد على منهجية للمسار التدريبى بحيث يكون المتدرب قد تلقى الدورات الحقلية بداية من التعرف على مراحل التنقية والمعالجة ثم الانتقال إلى الدورات المعملية داخل معمله طبقا للإطار الزمنى المحدد للمدد البينية لكل درجة وظيفية.

ولقد اعتمدنا على وضع معايير لكل مرحلة في إعداد هذا الاصدار وكان من أهم هذه المعايير:

- المشاركة الفعالة للخبرات والكفاءات التدريبية بالشركات التابعة في وضع المناهج بما يناسب عموم الكيميائيين على مستوى الجمهورية.
- عقد ورشة عمل متخصصة لكل مادة تدريبية يشارك بها جميع المدربين ذوى <u>التخصص والخبرات</u> سواء من المعمل المرجعى أو معامل الشركات فضلا عن أن يكون المدرب قد قام بتدريس هذه المادة مرات عديدة.
- 3. استخدام وسيلة اتصال غير تزامنى بين جميع المدربين المعتمدين لكل مادة على حدة من خلال انشاء جروب على الفيس بوك لكل مادة على حده (مذكور فى دليل المدرب).
- وضم حقيبة تدريبية كاملة لكل برنامج معدة طبقا لأحدث النظم والمعايير العالمية تحتوى على (دليل المتدرب- شرائح العرض- ملحقات مقرؤة ومرئية- دليل المدرب-بنك الأسئلة).
- 5. بناء المحتوى لكل برنامج تدريبى طبقاً لأحدث المراجع العالمية ومن أهمها كتاب الطرق القياسية لتحليل مياه الشرب والصرف الصحى (الإصدار رقم 23) وبما يتوافق مع متطلبات آخر إصدارات الايزو(1702)، مع مراعاة التحديثات الخاصة بالتشريعات والقوانين ذات الصلة.

وجدير بالذكر أن الإصدار الثانى من البرامج التدريبية اعتمد فى تصميمه على عرض مبسط للمعلومات قدر الامكان طبقاً للأسس العلمية وطبقاً للجدارات المعتمدة على تحديد أهداف واضحة وصريحة لتدريب المتدربين، وتشتق تلك الجدارات من الفهم الواضح لدور المتدرب طبقا لبطاقة الوصف الوظيفى، وتتضمن معارف ومهارات وسلوك. مما يضمن إكساب المتدرب مهارات سلوكية بالإضافة إلى المواد التخصصية. كما تم تصميم العديد من ورش العمل على أسماس تسميل و تسريع عمليتي التعلم و كسب المهارات بما يسمح بتعظيم الفائدة من العملية التدريبية.

كذلك تم استخدام أساليب التدريب الحديثة والاعتماد على التدريب التفاعلى والتركيز على الجوانب التطبيقية في استخدام الوسائل والأساليب المختلفة ، كما تم استخدام الطرق الحديثة للتعليم التفاعلي والغير تزامني كمصادر مساندة للتعلم من خلال انشاء جروب على الفيس بوك للمدربين المعتمدين (HCWW Trainers).

وفى الختام نرجوا من الله أن يتقبل منا هذا العمل كما نأمل أن يكون هذا العمل علما نافعا للعاملين بقطاع المعامل بالشـركة القابضـة والشـركات التابعة لما يشـمله من معلومات فنية قيمة وأن يفيد العاملين الجدد بها ليصبحوا قادرين على تنفيذ مهامهم الوظيفية بالشكل الأمثل .

والله ولى التوفيق.

1. INTRODUCTION

1. Definitions

- Water Microbiology: Study of small living microorganisms in water.
- **Drinking Water**: Regarding to World Health Organization (WHO) should be "Suitable for human consumption and for all usual domestic purposes".
- Drinking Water (from microbiological perspective): Should be free of any microorganism that may cause health risk (pathogen).
- **Pathogen**: Is an organism with the potential to cause disease.

2. Why Water microbiology is important?

- It helps us to:
- Identify harmful microorganisms in water.
- Realize the diseases cause.
- Measure effectiveness of water treatment.
- Monitor water safety.
- Confirm that water will not cause human disease.

3. History of Water Microbiology

- 1. 1852: British physician, Dr. John Snow (1813 1858), traced Cholera epidemic in London.
- 1. He noticed that cholera cases were isolated from region that served by specific company, while other areas did not affected.
- 2. After removal of the pump, no other cholera cases were found.
- 3. He concluded that "water can cause human health disease".

2. 1880, French chemist Louis Pasteur (1822-1895)

1. He defined the "germ theory" of disease, stating that "infectious disease wasn't spontaneous but a result of the presence and/or growth of harmful microorganisms".

3. 1877, German physician Robert Koch (1843-1910)

1. isolated Bacillus anthracis as the causative agent Anthrax.

4. 1883, Koch published an article entitled: "About Detection Methods for Microorganisms in Water"

- 1. In that historic paper he marked the introduction of the application of microbial indicators.
- 2. He described for the first time the methodology for bacteria measurement in water.

5. 1884, Koch identified *Vibrio cholerae* as the causative agent of cholera after his research with French team in governmental hospital in Alexandria, Egypt.

1. He is considered one of the founders of microbiology science.

6. 1850, Dr. Snow used Chlorine as water disinfection agent in England.

7. 1918, The US Department of Treasury called all world water companies to disinfect water with chlorine to kill harmful microorganisms.

1. Bacteriological Indicators

1. Concepts:

• **Microorganism** or group of microorganisms that its presence in water indicate the possibility of pathogens contamination,

التحاليل البكتيريولوجية

- Its absence indicates no pathogenic contamination. No single indicator organism can universally be used for all purposes of water quality surveillance.
- The challenge is to select the appropriate indicator, or combination of indicators, for each particular purpose of water quality assessment.

2. Why indicators and not direct detection of pathogens?

• We use indicators because:

- **Pathogens** are in small numbers.
- Takes long times to detect (72 hours for Typhoid).
- Cannot predict likelihood of infectious dose from average concentration.
- Many pathogens present, and each one requires different method.

3. Characteristics of Microbial Indicators:

- Present when there is a risk of contamination by pathogens;
- Should not multiply in environmental conditions under which pathogens cannot multiply;
- Should correlate with the degree of faecal contamination;
- Easy to enumerate and identify by simple methods;
- Specific for faecal or sewage pollution;
- At least as resistant as pathogens to conditions in natural water environments, and water treatment processes;
- Should have stable characteristics and give consistent reactions in these analyses;
- Present in the same or higher numbers than pathogens;
- Non-pathogenic;

4. Examples:

- Water Treatment and Quality:
- Heterotrophic Plate Count.
- Total Coliform
 - Fecal Pollution:
- Thermotolerant Coliforms (Fecal Coliform).
- Escherichia coli.
- Intestinal Enterococci.

7

Indicator	Characteristic	Source and Occurrence
Total coliforms	include a wide range of aerobic and facultatively anaerobic. Gram-negative, non-spore-forming bacilli capable of growing in the presence of high bile salts concentrations with the fermentation of lactose (by me β -galactosidase enzyme) and production of acid or aldehyde within 24 hours at 35–37 °C.	occur in both sewage and natural waters. Some of these bacteria are excreted in the faeces of humans and animals, but many coliforms are heterotrophic and able to multiply in water and soil environments
Thermotolerant Coliform	Subset of Total Coliform that that can ferment lactose at higher temperatures (44.5°C).	
E. coli	Rod-shaped; ferment lactose and produce gas at 44.5 C.	Occurs in high numbers in human and animal faeces, sewage and water subject to recent faecal pollution.
Intestinal Enterococci	Gram-positive cocci, relatively tolerant of sodium chloride and alkaline pH levels. Facultatively anaerobic and occur singly, in pairs or as short chains.	excreted in the faeces of humans and other warm-blooded animals. Some have also been detected in soil in the absence of faecal contamination, present in sewage and water environments polluted by sewage or wastes from humans and animals.



Relation between microbial indicators

8. Methods of Analysis

1. Estimation of Heterotrophic Bacteria by Pour Plate Technique

1. Scope:

- This method describes a detailed procedure for estimation of live culturable heterotrophic bacteria by pour plate method in water samples in 48 h at 35.0 °C and 7 days at 22.0 °C on the basis of culturing in non-selective high and low nutrient agar medium respectively.
- This procedure can be applied for different types of water e.g. raw water, ground water, network distribution system (drinking water), water from different treatment process...etc.

Tip: As pour plate can accommodate large sample volume, it is better to be applied in analysis of ground and treated water.

2. Principle:

- Volume of water to be tested is mixed with Plate Count Agar/ R2A media in sterile petri dishes, incubated at 35±0.5°C for 48h and/or 22±1.0°C for 7 days respectively.
- The advantage of R2A medium is due to low nutrient formulation that stimulates the growth of stressed and chlorine-tolerant bacteria. The final count also depends on interaction among the developing colonies.
- The test has little value as an indicator of pathogen presence but can be useful in operational monitoring as a treatment and disinfectant indicator.
- In addition, HPC measurement can be used as indicator of stagnation, tuberculation, residual disinfectant concentration, assessing the cleanliness and integrity of distribution systems, presence of biofilms and availability of nutrients for bacterial growth.
- It should be noted that heterotrophic plate count results are not an indicator of water safety and, as such, should not be used as an indicator of potential adverse human health effects. It should also be noted that the results obtained using heterotrophic plate count test are not an accurate assessment of total heterotrophic concentrations but, instead, are indications of culturable organisms present.
- Pour plate method is a simple technique that can accommodate volumes of sample or diluted sample ranging from 0.1 to 2.0 ml. The colonies produced are relatively

small and compact, showing fewer tendencies to encroach on each other than those produced by surface growth. On the other hand, submerged colonies often are slower growing and are difficult to transfer.

3. Definition:

- Heterotrophs are those microorganisms that use organic compounds for most or all of their carbon requirements.
- Heterotrophic Plate Count (HPC): Is a procedure for estimating the number of live cultivable heterotrophic bacteria in water and measuring changes during water treatment and distribution or in swimming pools.
- Note: Various synonyms are frequently used instead of "heterotrophic count". These include colony count, viable count, standard plate count and culturable micro-organisms.

4. Environmental Conditions:

• All sample analysis steps will be carried out under aseptic conditions (Laminar Air Flow "LAF" is preferred)

5. Interference:

• Significant heat shock to bacteria from the transient exposure of the sample to 45 to 46 °C agar may occur, to minimize, use a thermostatically controlled water bath for tempering the agar.

6. Equipment:

- Water bath: 45°C to 46°C.
- Sterile Petri dishes (preferably plastic) about 90 x 15 mm, with tight-fitting lids.
- Automatic pipettes and associated sterile tips capable of delivering 1 ml.
- Automatic pipettes and associated sterile tips capable of delivering 100 µl.
- Buffered Dilution water screw cap bottles.
- Incubator: $35^{\circ}C \pm 0.5^{\circ}C$.
- Incubator: $22^{\circ}C \pm 1^{\circ}C$.
- Colonies counter aid.

7. Chemicals and Reagents:

- Plate Count Agar.
- R2A medium
- Phosphate buffered dilution water.

8. **Precautions:**

- Follow the normal safety procedures required in a microbiology laboratory.
- Mouth-pipetting is prohibited.
- Do not re-sterilize plating medium, re-melting is allowed for only one time.
- Decontaminate all the used plates and materials at the end of the analyses.
- Do not prepare dilutions or pour plates in direct sunlight.
- Avoid spilling medium on outside of dish lid when pouring plates.
- Do not depend on the sense of touch to indicate proper medium temperature when pouring agar.
- Avoid prolonged exposure to unreasonably high temperatures during melting agar.
- Use sterile tips for initial and subsequent transfers from each different dilution.
- Replace the tip with a sterile one, if it becomes contaminated before transfers completed.
- Avoid contamination when removing sterile tips from the container by not dragging tip across exposed ends of tip in the tip container or across lips and necks of dilution bottles.
- Do not insert tips more than 2 to 3 cm below the surface of sample or dilution when removing the sample volume.
- Discard melted agar that contains precipitate.

9. **Procedure:**

1. Sample Handling:

- Initiate analysis as soon as possible after collection to minimize changes in bacterial populations, and do not exceed a holding time of 24 h.
- If the sample cannot be tested within 30 min after collection, maintain it at <10°C but do not allow it to freeze during transit.

2. Instrument Verification:

- Calibrate and verify balances daily using reference weights.
- standardize pH meter prior to use, using at least two standards points.
- Check temperature of incubators and refrigerators daily to ensure operation within stated limits of method.
- Check sterilization procedure and efficiency to ensure media sterility used in this method.

3. Media preparation:

• Prepare and handle media and reagents appropriately, carry out suitable performance test to ensure homogeneity, sterility, and suitability of each prepared batch.

4. Melting Agar:

- Melt sterile solid agar medium in boiling water bath.
- Maintain melted medium in a water bath between 44°C and 46°C until used, preferably for no longer than 3 h.

5. Sample Dilution and preparations:

- Label the bottom of plates with the sample identification number, analyst initials, date and sample volume to be analyzed.
- Prepare buffered dilution water according to Annex A.1.
- Select the dilution(s) so that the total number of colonies on a plate will be between 30 and 300. For example, where a heterotrophic plate counts high as 3000 is suspected, prepare plates with 10–2 dilution.
- Shake the sample container vigorously by rapidly making about 25 complete up-anddown (or back-and-forth) movements. Optionally, use a mechanical shaker to shake samples or dilutions for 15s.
- Hold micropipette inside neck of dilution bottle to transfer 10 ml of sample to 90 ml of sterile dilution water bottle (or 1 ml of sample to 09 ml of diluent), cap, and mix by vortex, 1 ml of this dilution is considered 10-1 of the original sample.
- Repeat the previous step to prepare further dilutions.

6. Selection of Sample Size:

• Pipette 1.0 ml, and 0.1 ml directly from tested water sample. In examining, known, turbid water, do not measure 0.1 ml inoculum of original sample, but prepare an appropriate dilution.

7. Plating:

- Limit the number of samples to be plated in any one series so that no more than 20 minutes (preferable 10 minutes) elapse between dilution of the first sample and pouring of the last plate in the series.
- Prepare two sets of plates for each sample; PCA set to be incubated at 35.0 °C, and R2A set to be incubated at 22.0 °C.
- Shake the sample container vigorously.



- Lift cover of petri dish just high enough to insert micropipette.
- Pipet selected sample volume by holding micropipette at an angle of about 45°.
- Remove micropipette without touching it to dish.



- Pour at least 10 to 12 ml liquefied medium maintained at 44 to 46°C into each dish by gently lifting cover just high enough to pour.
- Mix melted medium carefully with test portions in petri dish.
- Rotate the dish gently containing the mixture first in one direction and then in the opposite direction.
- Let plates solidify on a level surface.

8. Incubation:

- Invert plate count agar plates and incubate plates at 35±0.5°C for 48hr
- Invert R2A plates and incubate for 22±1.0°C for 7 days.
- Counting:
- Count colonies using dark-field colony counter.
- If counting must be delayed temporarily, store plates at 5 to 10°C for no more than 24 h.



Quebec counter

- Count all colonies on selected plates promptly after incubation.
- Count plates that have 30 to 300 colonies.
- If there is no plate with 30 to 300 colonies, and one or more plates have more than 300 colonies, use the plate(s) having a count nearest 300 colonies.
- When fewer than 10 colonies/cm2 are present, count colonies in 19 squares (of the colony counter) having representative colony distribution. Multiply sum of the number of colonies by 3 to compute estimated colonies per plate when the plate area is 57 cm2. if plate area is 65 cm2(Glass) multiply sum of the number of colonies in 13 square centimeters by 5.
- When more than 10 colonies/cm2 are present, count four representative squares, take average count per square centimeter, and multiply by 57 to compute estimated colonies per plate when the plate area is 57 cm2 or multiply by 65 when the plate area is 65 Cm2.
- When the total number of colonies is less than 30, disregard the rule above and record result observed.
- Count as individual colonies similar-appearing colonies growing in close proximity but not touching, provided that the distance between them is at least equal to the diameter of the smallest colony.



• Count impinging colonies that differ in appearance, such as morphology or color, as individual colonies.



• If spreading colonies (spreaders) are encountered on the plate(s) selected, count colonies on representative portions only when colonies are well distributed in spreader-free areas and the area covered by the spreader(s) does not exceed one-half the plate area.

• Count each of the following types of spreading colonies as one: a chain of colonies that appears to be caused by disintegration of a bacterial clump as agar and sample were mixed; a spreader that develops as a film of growth between the agar and bottom of Petri dish; and a colony that forms in a film of water at the edge or over the agar surface.



9. Calculations:

• Compute bacterial count per milliliter by the following equation:

$$CFU/mL = \frac{\text{colonies counted}}{\text{actual volume of sample in dish, mL}}$$

10. Quality Control:

1. Analytical Quality Control Batch:

- Prepare duplicate plates for each volume of sample or dilution examined; concentration value for that volume will be obtained by getting the arithmetic mean for both plates.
- Plate(s) with non-selective agar medium, labeled as "Environmental Control", will be opened along the samples analysis run time to ensure adequate environmental conditions.
- Procedural Blank: Inoculate 1 ml of sterile phosphate buffered dilution water, using separate sterile micropipette tip, after inoculation of samples series; incubate under

the same conditions as a sample. Absence of growth indicates absence of cross contamination.

11. Reporting:

- Report results as the number of colony forming units (CFU) per ml.
- If sample was delayed more than permissible time, write "Delayed" in report comment.
- If plates from all dilutions of any sample have no colonies, report the count as less than one (< 1) divided by the corresponding largest sample volume used. For example, if no colonies develop from the 0.01-ml sample volume, report the count as less than 100 (< 100) CFU/ml
- When colonies on duplicate plates and/or consecutive dilutions are counted and results are averaged before being reported, round off counts to two significant figures only when converting to colony-forming units.
- Raise the second digit to the next higher number when the third digit from the left is 5, 6, 7, 8, or 9; use zeros for each successive digit toward the right from the second digit. For example, report a count of 142 as 140 and a count of 155 as 160, but report a count of 35 as 35.
- If the number of colonies per plate far exceeds 300, do not reports result as "too numerous to count" (TNTC) otherwise use one of the following rules:
- When bacterial counts on crowded plates are greater than 100 colonies/cm2, report result as greater than (>) 5700 divided by the smallest sample volume plated for plastic plates. Report as estimated colony-forming units per milliliter.
- If plates have excessive spreader growth, report as "spreaders" (Spr)
- When plates are uncountable because of missed dilution, accidental dropping, and contamination, or the control plates indicate that the medium or other material or lab ware was contaminated, report as "laboratory accident" (LA).

12. References:

- Eugene W. Rice, Rodger B. Baired, Andrew D. Eaton, Lenore S. Clesceri eds.2017 Standard Methods for the Examination of Water and Wastewater. Method 9215B#, 23rd edition. American Public Health Association, American Water Works Association, and Water Environment Federation
- World Health Organization. 2011, Guidelines for Drinking- Water Quality, Fourth Edition. Geneva,

2. Estimation of Heterotrophic Bacteria by Spread Plate Method

1. Scope:

- This method describes a detailed procedure for estimation of live culturable heterotrophic bacteria by spread plate method in water samples in 48 h at 35.0 °C and 7 days at 22.0 °C on the basis of culturing in non-selective high and low nutrient agar medium respectively.
- This procedure can be applied for different types of water e.g. raw water, ground water, network distribution system (drinking water), water from different treatment process...etc.

Tip: As no regulation present for raw water, it is better to apply spread plate technique for better recovery of bacterial population

2. Principle:

- Volume of water to be tested is spread over medium in sterile petri dishes, and incubated at 35°C for 48h and/or 22°C for 7 days respectively.
- The advantage of R2A medium is due to low nutrient formulation that stimulates the growth of stressed and chlorine-tolerant. The final count also depends on interaction among the developing colonies.
- The test has little value as an indicator of pathogen presence but can be useful in operational monitoring as a treatment and disinfectant indicator. In addition, HPC measurement can be used as indicator of stagnation, tuberculation, residual disinfectant concentration, assessing the cleanliness and integrity of distribution systems, presence of biofilms and availability of nutrients for bacterial growth.
- It should be noted that heterotrophic plate count results are not an indicator of water safety and, as such, should not be used as an indicator of potential adverse human health effects. It should also be noted that the results obtained using heterotrophic plate count test are not an accurate assessment of total heterotrophic concentrations but, instead, are indications of culturable organisms present.
- Spread plate method has the advantage of using solidified agar, eliminating the possibility of heat shock. The resultant colonies can be easily transferred, and their colony morphology can be distinguished. However, this method is limited by the small volume of sample or diluted sample that can be absorbed by the agar: 0.1 to 0.5 ml, depending on the degree to which the pre-poured plates have been dried.

3. Definitions:

- Heterotrophs are those microorganisms that use organic compounds for most or all of their carbon requirements.
- Heterotrophic plate count (HPC): is a microbial method that uses colony formation on culture media to approximate the levels of heterotrophic flora.
- Colony-Forming Units (CFU): Colonies may arise from pairs, chains, clusters, or single cells.
- Note: Various synonyms are frequently used instead of "heterotrophic count". These include colony count, viable count, standard plate count and culturable micro-organisms etc.

4. Interference:

• Spreading colonies, such as a film of growth between the agar and bottom of petri dish or a film of water at the edge or over the agar surface, are mainly develop because of an accumulation of moisture at the point from which the spreader originates. If those spreading colonies cover more than half the plate, it will interfere with obtaining reliable plate count.

5. Equipment:

- Pre-sterilized plastic Petri dishes about 90 x 15 mm, with tight-fitting lids.
- Automatic pipettes and associated sterile tips capable of delivering 100 µl.
- Automatic pipettes and associated sterile tips capable of delivering 1 ml
- Pre-sterilized plastic spreader.
- Buffered dilution screw-capped bottles.
- Incubator: $35^{\circ}C \pm 0.5^{\circ}C$.
- Colonies counter aid.
- Indelible ink marker for labeling plates.

6. Chemicals and Reagents:

- Plate Count Agar.
- R2A medium
- Phosphate buffered dilution water.

7. Precautions:

• Follow the normal safety procedures required in a microbiology laboratory.

- Mouth-pipetting is prohibited.
- Do not re-sterilize plating medium, re-melting is allowed for only one time..
- Decontaminate all the used plates and materials at the end of the analyses.
- Do not prepare dilutions or pour plates in direct sunlight.
- Avoid prolonged exposure to unreasonably high temperatures during melting agar.
- Use sterile tips for initial and subsequent transfers from each different dilution.
- Replace the tip with a sterile one, if it becomes contaminated before transfers completed.
- Avoid contamination when removing sterile tips from the container by not dragging tip across exposed ends of tip in the tip container or across lips and necks of dilution bottles.
- Do not insert tips more than 2 to 3 cm below the surface of sample or dilution when removing the sample volume.
- Discard melted agar that contains precipitate.
- Pre-dry the dish of excess moisture before use to prevent water film formation.
- Allow the inoculum to soak into the agar to prevent colony spreading over agar layer.

8. Procedure:

1. Sample Handling:

- Initiate analysis as soon as possible after collection to minimize changes in bacterial populations, and do not exceed a holding time of 24 h.
- If the sample cannot be tested within 30 min after collection, maintain it at <10°C but do not allow it to freeze during transit.

2. Media and solutions preparation:

- Prepare and handle media and reagents appropriately, carry out suitable performance test to ensure homogeneity, sterility, and suitability of each prepared batch.
- Phosphate buffered dilution water: See annex A.1.
- If plates are previously prepared and stored in the refrigerator, allow them to warm at room temperature. The crystals that form on agar plates after refrigeration will disappear as the plates warm up.
- Use plates within two weeks after preparation.

3. Preparation of Plates:

- Pour 15 ml of medium into sterile 90×15 petri dishes.
- Let agar solidify.

- Use pre-dried plates immediately after drying or store for up to 2 weeks refrigerated at 2 -8 °C.
- If pre-drying and using plates the same day, pour 25 ml agar into petri dish and dry in laminar air flow at room temperature (24°C-26 °C) with the lid off to obtain the desired 2- to 3-g weight loss.

4. Sample Dilution and preparations:

- Label the bottom of plates with the sample identification number, analyst initials, and sample volume to be analyzed.
- Prepare buffered dilution water as directed in Annex 01.
- Select the dilution(s) so that the total number of colonies on a plate will be between 30 and 300. For example, where a heterotrophic plate counts high as 3000 is suspected, prepare plates with 10–2 dilution.
- Use decimal dilutions in preparing sample volumes of less than 0.1 ml.
- Shake the sample container vigorously by rapidly making about 25 complete up-anddown (or back-and-forth) movements. Optionally, use a mechanical shaker to shake samples or dilutions for 15 s.
- Hold micropipette inside neck of sample bottle to transfer 1.0 ml of sample to 9 ml of sterile dilution water tube (or 10 ml to 90 ml of diluent), cap, and mix by vortex.
- $100 \ \mu l$ of this dilution is considered 10-2.
- Repeat the previous step to prepare further dilutions.

5. Selection of Sample Size:

• For raw water, Pipette (500 μ l and 100 μ l) from 10-1 dilution, or as appropriate.

6. Inoculation:

- Allow the dish to equilibrate to room temperature when a prepared petri dish is used.
- Pipet selected volume onto surface of pre-dried agar plates.
- Distribute inoculum over surface of the medium, using a sterile spreader, by manually rotating.
- Let inoculum to be absorbed completely onto the medium before incubation.

7. Incubation:

- Invert PCA plates and incubate at 35±0.5°C for 48 hr.
- Invert R2A plates (if any) and incubate at 22.0±1.0°C for 7 days.

8. Counting

• Count colonies using dark-field colony counter.

- If counting must be delayed temporarily, store plates at 5 to 10°C for no more than 24 h.
- Count all colonies on selected plates promptly after incubation.
- Count plates that have 30 to 300 colonies.
- If there is no plate with 30 to 300 colonies, and one or more plates have more than 300 colonies, use the plate(s) having a count nearest 300 colonies.
- When the total number of colonies is less than 30, disregard the rule above and record result observed.
- Count as individual colonies similar-appearing colonies growing in close proximity but not touching, provided that the distance between them is at least equal to the diameter of the smallest colony.
- Count impinging colonies that differ in appearance, such as morphology or color, as individual colonies.
- If spreading colonies (spreaders) are encountered on the plate(s) selected, count colonies on representative portions only when colonies are well distributed in spreader-free areas and the area covered by the spreader(s) does not exceed one-half the plate area.
- Count each of the following types of spreading colonies as one: a chain of colonies that appears to be caused by disintegration of a bacterial clump as agar and sample were mixed; a spreader that develops as a film of growth between the agar and bottom of Petri dish; and a colony that forms in a film of water at the edge or over the agar surface

9. Calculations:

• Compute bacterial count per milliliter by the following equation:

 $CFU/mL = \frac{\text{colonies counted}}{\text{actual volume of sample in dish, mL}}$

9. Quality Control:

• Prepare duplicate plates for each volume of sample or dilution examined; concentration value for that volume will be obtained by getting the arithmetic mean for both plates.

- Plate(s) with non-selective agar medium, labeled as "Environmental Control", will be opened along the samples analysis run time to ensure adequate environmental conditions.
- Procedural Blank: Inoculate 1 ml of sterile phosphate buffered dilution water by pour plate technique, using separate sterile micropipette tip, after inoculation of samples series; incubate under the same conditions as a sample. Absence of growth indicates absence of cross-contamination.

10. Reporting:

- Report results as the number of colony forming units (CFU) per ml.
- If sample was delayed more than permissible time, write "Delayed" in report comment.
- If plates from all dilutions of any sample have no colonies, report the count as less than one (< 1) divided by the corresponding largest sample volume used. For example, if no colonies develop from the 0.01-ml sample volume, report the count as less than 100 (< 100) CFU/ml
- When colonies on duplicate plates and/or consecutive dilutions are counted and results are averaged before being reported, round off counts to two significant figures only when converting to colony-forming units.
- Raise the second digit to the next higher number when the third digit from the left is 5, 6, 7, 8, or 9; use zeros for each successive digit toward the right from the second digit. For example, report a count of 142 as 140 and a count of 155 as 160, but report a count of 35 as 35.
- If the number of colonies per plate far exceeds 300, do not report result as "too numerous to count" (TNTC) otherwise use one of the following rules:
- When fewer than 10 colonies/cm² are present, count colonies in 19 squares (of the colony counter) having representative colony distribution. Multiply sum of the number of colonies by 3 to compute estimated colonies per plate when the plate area is 57 cm2.
- When more than 10 colonies/cm² are present, count four representative squares, take average count per square centimeter, and multiply by 57 to compute estimated colonies per plate when the plate area is 57 cm².

- When bacterial counts on crowded plates are greater than 100 colonies/cm², report result as greater than (>) 5700 divided by the smallest sample volume plated for plastic plates. Report as estimated count /ml.
- If plates have excessive spreader growth, report as "spreaders" (Spr)
- When plates are uncountable because of missed dilution, accidental dropping, and contamination, or the control plates indicate that the medium or other material or lab ware was contaminated, report as "laboratory accident" (LA).

11. References:

 Eugene W. Rice, Rodger B. Baired, Andrew D. Eaton, Lenore S. Clesceri eds.2017 Standard Methods for the Examination of Water and Wastewater. Method 9215C#, 23rd edition. American Public Health Association, American Water Works Association, and Water Environment Federation

3. Estimation of Coliform Group by Multiple Tube Fermentation Technique (For Drinking Water)

1. Scope:

- The procedure below is for estimation of Coliforms by Multiple-Tube Fermentation technique (MTF), also called Most Probable Number (MPN) procedure, in 96 hours, or less, in water samples on the basis of the production of gas and acid from fermentation of lactose.
- This procedure can be applied for potable water.

2. **Principle:**

- Volumes of water sample to be tested are added to tubes, or bottles, containing the presumptive media with inverted vials (Durham Tube) or pH indicator and incubated at 35°C for 48 hours. The selectivity of media is due to sodium lauryl sulfate that acting as inhibitor of bacteria other than coliforms. After incubation, the tubes, or bottles, are examined for growth, gas, and/or acidic reaction (shades of yellow color).
- An additional confirmatory test is required to confirm the result.
- Total coliform bacteria (excluding *E. coli*) occur in both sewage and natural waters. Some of these bacteria are excreted in the faeces of humans and animals, but many coliforms are heterotrophic and able to multiply in water and soil environments.
- Total coliforms include organisms that can survive and grow in water, particularly in the presence of biofilms. They were, traditionally, regarded as belonging to the genera *Escherichia, Citrobacter, Klebsiella and Enterobacter*, but the group is more heterogeneous and includes a wider range of genera, such as Serratia and Hafnia. Hence, they are not useful as an index of faecal pathogens, but they can be used as an indicator of treatment effectiveness and to assess the cleanliness and integrity of distribution systems and the potential presence of biofilms. However, there are better indicators for these purposes.
- Thermotolerant coliform (Fecal coliform) include species that may inhabit the intestines of warm-blooded animals and human. In most waters, the predominant genus is Escherichia, but some types of *Citrobacter, Klebsiella and Enterobacter* are also thermotolerant. They are usually found in sewage and water recently subjected to fecal pollution.

• Populations of thermotolerant coliforms are composed predominantly of *E. coli*; as a result, this group is regarded as a less reliable but acceptable index of faecal pollution, their presence in drinking water indicates the possible presence of pathogens.

		β-Galactosidase
Escherichia	Escherichia	Escherichia
Klebsiella	Klebsiella	Klebsiella
Enterobacter	Enterobacter	Enterobacter
Citrobacter	Citrobacter	Citrobacter
	Yersinia	Yersinia
	Serratia	Serratia
	Hafnia	Hafnia
	Pantoea	Pantoea
	Kluyvera	Kluyvera
		Cedecea
		Ewingella
		Moellerella
		Leclercia
		Rahnella
		Yokenella

bold type = coliforms which can be present in the environment as well as in human faeces.

bold and underline = coliforms which are considered to be primarily environmental.

Total Coliform species according to testing method

• Results are reported in terms of the Most Probable Number (MPN) of organisms present. This number, based on certain probability formulas, is an estimate of the mean density of total coliforms in the sample.

3. Definitions:

- Total coliform bacteria in this procedure are those facultative anaerobic, gramnegative, non-spore-forming, rod-shaped bacteria that ferment lactose with gas and acid formation within 48 hours at 35°C.
- Thermotolerant coliform (Fecal coliform) bacteria in this procedure are those facultative anaerobic, gram-negative, non-spore-forming, rod-shaped bacteria that ferment lactose with gas and acid formation within 48 hours at 44.5°C.

4. Environmental Conditions:

• All sample analysis steps will be carried out under aseptic environmental conditions.

5. **Precautions:**

- Follow the normal safety procedures required in microbiology laboratory.
- Mouth-pipetting is prohibited.
- Use a separate sterile tip for each sample.
- Replace the tip with a sterile one if becomes contaminated before transfers are completed.
- Do not insert tips more than 2.5 cm below the surface of sample when removing the sample volume.
- Avoid contamination when removing sterile tips from the container by not dragging tip across lips and necks of dilution bottles.
- Avoid picking up any membrane or scum on the needle while transferring by inclining the fermentation tube.
- Insert end of loop or needle into the liquid in the tube to a depth of approximately 0.5 cm.
- Decontaminate all the used tubes, plates, and materials at the end of the analyses.

6. Interference:

- Since the MPN indexes are based on a Poisson distribution, if the sample is not adequately mixed to ensure equal bacterial cell distribution before portions are removed, the MPN value will be a misrepresentation of the bacterial density.
- High densities of non-coliform bacteria and the inhibitory nature of some MTF media, e.g. Lauryl Tryptose broth, may prevent gas formation. It's recommended to treat all tubes with turbidity, indicating growth, regardless of gas production, as presumptive total coliform-positive tubes.

7. Equipment:

- 16×100 mm tubes, borosilicate glass, with heat-resistant plastic caps.
- 16×150 mm tubes, borosilicate glass, with heat-resistant plastic caps.
- Durham tubes, borosilicate glass.
- Test tubes racks to hold sterile culture tubes.
- Automatic pipettes and associated sterile long tips capable of delivering the selected sample volume.
- Incubator: $35^{\circ}C \pm 0.5^{\circ}C$.
- Water bath Incubator: $44.5^{\circ}C \pm 0.2^{\circ}C$.

- Sterile disposable applicator stick.
- Sterile glass/ plastic petri dishes 90 mm.
- Sterile glass/ plastic petri dishes 50 mm.
- Indelible ink marker for labeling plates.

8. Chemicals and Reagents:

- Lauryl Tryptose broth.
- Brilliant Green broth.
- Ec broth
- Phosphate buffered dilution water.
- Nutrient agar.
- Gram Stain Reagent set.
- M-Endo LES agar.
- Bromcresol purple.

9. Procedure:

1. Sample Handling:

- Analyze samples on day of receipt whenever possible. The time between sample collection and the placement of samples in the incubator must not exceed 30 hours.
- If arrival is too late for processing on same day, refrigerate overnight as long as holding time conditions still be met i.e. preservation temperature <10.0°C and time did not exceeded 30 h.

2. **Presumptive Phase:**

• Selection of sample size:

• Inoculate a row of tubes with five 20-mL portions, ten 10-mL portions, or one 100 ml portion of water sample (for Presence-Absence test), based on the probable coliform density.

• Media Preparation:

- Store dehydrated media (powders) in tightly closed bottles in the dark at less than 30°C in an atmosphere of low humidity, e.g., desiccator. Do not use them if they are discolored or caked, or if the character of a free-flowing powder is lost.
- Prepare and handle media and reagents appropriately, and carry out suitable performance test to ensure homogeneity, sterility, and suitability for each prepared batch.

- Determine pH of medium after sterilization. The required final pH is given in the directions for preparing each medium.
- Prepare in accordance with following Table:

Inoculum mL	Amount of Medium in Tube <i>mL</i>	Volume of Medium + Inoculum mL	Dehydrated Laury Tryptose Broth Required <u>g/L</u>
1	10 or more	11 or more	35.6
10	10	20	71.2
10	20	30	53.4
20	10	30	106.8
100	50	150	106.8
100	35	135	137.1
100	20	120	213.6

• Lauryl Tryptose broth:

- Prepare and sterilize medium according to manufacturer instructions.
- Dispense the medium as follow:
- Dispense 50 ml of medium (strength 3x) in case of inoculating one 100 ml sample portion.
- Dispense 10 ml of medium (strength 3x) in case of inoculating five 20-mLsample portions.
- Dispense 10 ml of medium (strength 2x) in case of inoculating ten 10-mLsample portions.
- Place an inverted Durham tube inside the fermentation tubes before sterilization. Alternatively add 0.01 g/L bromcresol purple to medium to determine acid production.
- Close tubes (bottles) with heat-resistant plastic caps.
- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.
- Use tubes (bottles) within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes (bottles).
- Inoculation:
- Use lauryl tryptose broth for presumptive phase.
- Arrange fermentation tubes in rows of five or ten tubes each in a test tube rack.
- Shake sample and dilutions vigorously about 25 times.

- Inoculate each tube in the set with sample volume.
- Mix test portions in the medium by gentle agitation.
- Incubation:
- Incubate inoculated tubes (bottles) at $35 \pm 0.5^{\circ}$ C.
- After 24 ± 2 h swirl each tube (bottle) gently and examine it for growth, gas, and acidic reaction (shades of yellow color)
- If no gas or acidic reaction is evident, re-incubate and reexamine at the end of 48±3
 h.
- Record presence or absence of growth, gas, and acid production.
- If the inner vial is omitted, growth with acidity signifies a positive presumptive reaction.
- Interpretation:
- Production of an acidic reaction and/or gas in the tubes or bottles within 48 ± 3 h constitutes a positive presumptive reaction.
- Submit tubes or bottles with a positive presumptive reaction to the confirmed phase.
- The absence of acidic reaction and/or gas formation at the end of 48 ± 3 h of incubation constitutes a negative result.
- Submit drinking water samples demonstrating growth without a positive gas or acid reaction to the confirmed phase.

3. Confirmed Phase:

- Media Preparation:
- Brilliant green broth (For Total Coliform)
- EC broth (For Thermotolerant Coliform)
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml medium, before sterilization, in fermentation tubes with an inverted Durham tube.
- Close tubes with heat-resistant plastic caps.
- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes) after preparation.
- Inoculation:

- Arrange brilliant green and EC tubes rows, in separate racks, in a similar manner to positive presumptive tubes (bottles).
- Gently rotate positive presumptive tubes (bottle) to resuspend the organisms.
- Use a sterile loop 3.0 mm in diameter to transfer three loopfuls from positive presumptive tube to a brilliant green tube and EC tube for confirmation.
- Repeat for all other positive presumptive tubes.

• Incubation:

- Incubate Brilliant Green (Total Coliform) tubes at $35 \pm 0.5^{\circ}$ C for 48 ± 3 h.
- Incubate EC (Fecal Coliform) tubes at 44.5 ± 0.2 °C for 24 ± 2 h.

• Interpretation:

- Formation of gas in any amount, in Durham tube, of the brilliant green tubes, at any time, within 48 ± 3 h constitutes a positive confirmed result.
- Formation of gas in any amount, in Durham tube, of the EC tubes, at any time, within 24 ± 2 h constitutes a positive confirmed result.
- Calculate the MPN value of the number of positive brilliant green lactose bile tubes and EC tubes from MPN index see A.2.
- In case of inoculating one bottle with 100 ml sample portion, report result as present or absent.

4. Completed Phase (In case of Analyzing Total Coliform Only):

• Media Preparation:

- Use LES Endo agar plates and Nutrient agar slants for completed phase.
- m-Endo LES agar:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense in 5 to 7 mL quantities into lower section of 50 mm Petri dishes, if dishes of any other size are used, adjust quantity to give an equivalent depth of 4 to 5 mm.
- If plates are previously prepared and stored in the refrigerator, allow them to warm at room temperature. The crystals that form on agar plates after refrigeration will disappear as the plates warm up.
- Use plates within two weeks after preparation.

• Nutrient agar:

- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml in screw-capped tubes before sterilization.

- Place tubes in an inclined position immediately after sterilization, so that the agar will solidify with a sloped surface.
- Tighten screw caps after cooling and store in a protected, cool storage area.
- If plates or tubes are made ahead of time and stored in the refrigerator, remove them and allow warming at room temperature before use.

1. Procedure:

- Streak one LES Endo agar plate from each tube of brilliant green lactose bile broth showing gas, as soon as possible after the observation of gas.
- Change loop between second and third quadrants to improve colony isolation.
- Streak plates in a manner to insure presence of some separated discrete colonies.
- Incubate plates (inverted) at $35 \pm 0.5^{\circ}$ C for 24 ± 2 h.
- Pick from each plate one or more typical (pink to dark red with a green metallic surface sheen) well-isolated coliform colonies.
- If no typical colonies are present, pick two or more colonies considered most likely to consist of organisms of the coliform group.
- Transfer growth from each isolate to a single-strength lauryl tryptose broth fermentation tube and onto a nutrient agar slant.
- Incubate nutrient agar slants at 35 ± 0.5 °C for 24 ± 2 h.
- Incubate secondary broth tubes at 35 ± 0.5°C for 24 ± 2 h; if gas is not produced within 24 ± 2 h reincubate and examine again at 48 ± 3 h.
- Perform Gram staining technique from nutrient agar slant cultures corresponding to the secondary tubes that showed gas.
- Gram Stain technique:
- Prepare separate light emulsions of the test bacterial growth using drops of distilled water on the slide.
- Prepare separate light emulsions of the positive and negative control cultures on the same slide using drops of distilled water on the slide.
- Fix the slides by air drying.
- Stain for 1 min with ammonium oxalate-crystal violet solution, rinse slides in tap water and drain off excess.
- Apply Lugol's solution for 1 min, rinse stained slides in tap water.
- Decolorize for approximately 15 to 30 s with acetone alcohol by holding slide between the fingers and letting acetone alcohol flow across the stained smear until

the solvent flows colorlessly from the slide. Do not over-decolorize to not remove crystal violet from gram positive bacteria.

- Counterstain with safranin for 15 s, rinse with tap water.
- Blot dry with absorbent paper or air dry, and examine microscopically.
- Gram-positive organisms are blue (violet); gram-negative organisms are red.
- Results are only acceptable when controls given proper reactions.

• Interpretation:

• Formation of gas in the secondary tube of lauryl tryptose broth within 48 ± 3 h and demonstration of gram-negative, nonspore-forming, rod-shaped bacteria from the agar culture constitute a positive result for completed test.

1. Calculations:

- The MPN values, for variety of positive and negative tubes combinations, are given in appendix A.2.
- The samples volumes indicated in indexes A.2 illustrates MPN values for concentrations of Positive and Negative results when Five 20-ml or Ten 10-ml portions are used.
- If the sample portion volumes used are those found in the tables, report the value corresponding to the number of positive and negative results in the series as the MPN/100ml.

2. Quality Control:

- Apply the presumptive-confirmed phase of the multiple-tube procedure to all samples, if samples were analyzed for fecal coliform using MTF technique there will be no need to apply completed phase, otherwise proceed completed test to not less than 10% of all coliform-positive samples on a seasonal basis.
- Procedural Blank: perform, under the same conditions as sample, sterile phosphate buffered dilution water after inoculation of samples series. Absence of growth indicates sterility of inoculation series.

3. Reporting:

- Report coliform concentration as the Most Probable Number (MPN)/100 ml.
- Report coliform presence or absence as the "present" or "absent" (MPN)/100 ml.

كيميائي مياه – الدرجة الثالثة

4. References

- Eugene W. Rice, Rodger B. Baired, Andrew D. Eaton, Lenore S. Clesceri eds.2017 Standard Methods for the Examination of Water and Wastewater. Method 9221B#, 23rd edition. American Public Health Association, American Water Works Association, and Water Environment Federation
- World Health Organization. 2011, Guidelines for Drinking- Water Quality, Fourth Edition. Geneva
- Guidelines for Canadian Drinking water: Guideline Technical Document. 2006.
 Total Coliform. Health Canada. Ottawa. Ontario.

2. Estimation of Coliform Group by Multiple Tube Fermentation (For water other than Drinking Water)

1. Scope:

- Estimation of Total Coliforms (TC) by Multiple-Tube Fermentation technique (MTF), also called Most Probable Number (MPN) procedure, in 96 hours, or less, in water samples on the basis of the production of gas and acid from fermentation of lactose.
- This procedure can be applied for surface raw water.

2. Principle:

- Volumes/dilutions of water sample to be tested are added to tubes containing lauryl Tryptose broth as presumptive media (5 tubes per dilution) with inverted vials (Durham Tube), or pH indicator. The selectivity of media is due to sodium lauryl sulfate that acting as inhibitor of bacteria other than coliforms. After incubation, the tubes are examined for growth, gas, and/or acidic reaction (shades of yellow color) if inner vial is omitted, growth with acidity (yellow color) is signifies a positive presumptive reaction.
- An additional confirmatory test by brilliant green bile broth 2% is required to confirm the result. The selectivity is due to presence of both ox gall (bile) and brilliant green dye acting as inhibitor of gram positive and selected gram negative bacteria. Organisms, primarily coliform, which are resistant to those inhibitor and ferment lactose with gas formation, indicated by Durham tube, can replicate in this medium.
- Total coliform bacteria (excluding E. coli) occur in both sewage and natural waters. Some of these bacteria are excreted in the faeces of humans and animals, but many coliforms are heterotrophic and able to multiply in water and soil environments.
- Total coliforms include organisms that can survive and grow in water distribution system, particularly in the presence of biofilms. They were, traditionally, regarded as belonging to the



genera Escherichia, Citrobacter, Klebsiella and Enterobacter, but the group is more heterogeneous and includes a wider range of genera, such as Serratia and Hafnia. Hence, they are not useful as an index of faecal pathogens, but they can be used as an indicator of treatment effectiveness and to assess the cleanliness and integrity of distribution systems and the potential presence of biofilms. However, there are better indicators for these purposes.

- The presence of Total Coliform after disinfection indicates inadequate treatment, and presence in distribution system or stored water supplies can reveal regrowth and possible biofilm formation or contamination through ingress of foreign material, including soil or plants.
- The MTF procedure, in comparison with the membrane filter (MF) procedure, is more difficult to perform, takes longer to produce results, and lacks precision; the precision depends on the number of tubes used. The most satisfactory information will be obtained when the largest sample inoculum examined shows acid and/or gas in some or all of the tubes and the smallest sample inoculum shows no acid and/or gas in any or a majority of the tubes. However, the MTF procedure is still of value when conditions render the MF technique unusable e.g. with turbid, colored, or grossly contaminated water
- Results are reported in terms of the Most Probable Number (MPN/100ml) of organisms present. It can be estimated by the formula given below (section 11) or from MPN tables (Appendix 2) using the number of positive tubes in the multiple dilutions. This number, based on certain probability formulas, is an estimate of the mean density of total coliforms in the sample.

3. Definitions:

- Total coliform bacteria in this SOP are those facultative anaerobic, gram-negative, non-spore-forming, rod-shaped bacteria that ferment lactose with gas and acid formation within 48 hours at 35°C.
- Fecal coliform bacteria in this procedure are those facultative anaerobic, gramnegative, non-spore-forming, rod-shaped bacteria that ferment lactose with gas and acid formation within 48 hours at 44.5°C.

4. Environmental Conditions:

• All sample analysis steps will be carried under aseptic conditions.

5. **Precautions:**

- Follow the normal safety procedures required in microbiology laboratory.
- Mouth-pipetting is prohibited.

- Use a separate sterile tip for initial and subsequent transfers from each tube, and different dilution.
- Replace the tip with a sterile one if becomes contaminated before transfers are completed.
- Do not prepare dilutions in the direct of sunlight.
- Do not insert tips more than 2.5 cm below the surface of sample or dilution when removing the sample volume.
- Avoid contamination when removing sterile tips from the container by not dragging tip across lips and necks of dilution bottles.
- Avoid picking up any membrane or scum on the needle while transferring by inclining the fermentation tube.
- Insert end of loop or needle into the liquid in the tube to a depth of approximately 0.5 cm.
- Decontaminate all the used tubes, plates, and materials at the end of the analyses.

6. Interference:

- Since the MPN indexes are based on assumption of a Poisson distribution (random distribution) that's for if the sample is not adequately shaken before portions are removed or if bacterial cells clump, the MPN value will be an underestimate of actual bacterial density.
- High densities of non-coliform bacteria and the inhibitory nature of some MTF media, e.g. Lauryl Tryptose broth, may prevent gas formation. It's recommended to treat all tubes with turbidity, indicating growth, regardless of gas production, as presumptive total coliform-positive tubes in case if tubes form the next dilution gives true positive reaction.

7. Equipment:

- 16×100 mm tubes, borosilicate glass, with heat-resistant plastic caps.
- 16×150 mm tubes, borosilicate glass, with heat-resistant plastic caps.
- Durham tubes, borosilicate glass.
- Test tubes racks to hold sterile culture tubes.
- Automatic pipettes and associated sterile long tips capable of delivering 1 ml.
- Dilution screwed cap bottles.
- Incubator: $35^{\circ}C \pm 0.5^{\circ}C$.
- Circulating Water bath Incubator (preferably with gabled cover): $44.5^{\circ}C \pm 0.2^{\circ}C$.
- Sterile disposable applicator stick.
- Sterile plastic petri dishes 90 mm.
- Sterile plastic petri dishes 50 mm.
- Indelible ink marker for labeling plates.

8. Chemicals and Reagents:

- Lauryl Tryptose broth.
- Brilliant Green broth.
- Ec broth
- Phosphate buffered dilution water.
- Nutrient agar.
- Gram Stain Reagent set.
- M-Endo LES agar.
- Bromcresol purple.

9. Procedure:

1. Sample Handling:

- Analyze samples on day of receipt whenever possible. The time between sample collection and the placement of samples in the incubator must not exceed 30 hours.
- If arrival is too late for processing on same day, refrigerate overnight as long as holding time conditions still be met i.e. preservation temperature <10.0°C and time did not exceeded 30 h.

2. **Presumptive Phase:**

• Selection of sample size:

• Inoculate a five series of tubes with five dilutions of water sample (multiples and submultiples of 10 ml), based on the probable total coliform density.

• Media Preparation:

• Prepare and handle media and reagents appropriately, and carry out suitable performance test to ensure homogeneity, sterility, and suitability for each prepared batch.

• Phosphate buffered dilution water:

- See annex 01
- Lauryl Tryptose broth:

- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml medium, before sterilization, in fermentation tubes with an inverted Durham tube.
- Alternatively, and only in case of unavailability of Durham tubes, add 0.01 g/L bromcresol purple to medium to determine acid production.
- Close tubes with heat-resistant plastic caps.
- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes).
- Samples Dilutions:
- Mix the sample by vigorously shakes the bottle.
- Use a sterile tip to transfer 10 ml of sample to 90 ml, or 1 ml of sample to 9 ml of sterile dilution water bottle, cap, and mix. 1 ml of this dilution is considered 10-1 of the original sample.
- Repeat the previous step to prepare further dilutions.
- Inoculation:
- Use lauryl tryptose broth for presumptive phase.
- Arrange fermentation tubes in rows.
- Prepare five sets of tubes (e.g.10,1.0, 0.1, 0.01, 0.001 ml) using five tubes per set.
- Shake sample or dilutions vigorously.
- Inoculate each tube in the set with sample or dilution volume.
- Mix test portions in the medium by gentle agitation.
- Incubation:
- Incubate inoculated tubes at $35 \pm 0.5^{\circ}$ C.
- After 24 ± 2 h swirl each tube gently and examine it for growth, gas, and acidic reaction (shades of yellow color). If inner vial is omitted, growth with acidity (yellow color) is signifies a positive presumptive reaction
- If no gas or acidic reaction is evident, reincubate and reexamine at the end of 48 ± 3 h.
- Record presence or absence of growth, gas, and acid production.



• Interpretation:

- Production of an acidic reaction and/or gas in the tubes within 48 ± 3 h constitutes a positive presumptive reaction.
- Submit tubes with a positive presumptive reaction to the confirmed phase.
- The absence of acidic reaction and/or gas formation at the end of 48 ± 3 h of incubation constitutes a negative test.
- **3. Confirmed Phase:**
- Media Preparation:
- Brilliant green broth (For Total Coliform):
- EC broth (For Thermotolerant Coliform)
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml medium, before sterilization, in fermentation tubes with an inverted Durham tube.
- Close tubes with heat-resistant plastic caps.
- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.

- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes) after preparation.
- Inoculation:
- Submit tubes with a positive presumptive reaction to the confirmed phase within 24 ± 2 h of incubation to the confirmed phase.
- If additional presumptive tubes showed active fermentation or acidic reaction at the end of 48 ± 3 h incubation period, submit to the confirmed phase.
- Arrange brilliant green and EC tubes rows, in separate racks, in a similar manner to positive presumptive tubes (bottles).
- Gently rotate positive presumptive tubes (bottle) to resuspend the organisms.
- Use a sterile loop 3.0 mm in diameter to transfer three loopfuls from positive presumptive tube to a brilliant green tube and EC tube for confirmation.
- Repeat for all other positive presumptive tubes.
- Repeat for all other positive presumptive tubes.
- Incubation:
- Incubate the inoculated tubes at $35 \pm 0.5^{\circ}$ C for 48 ± 3 h
- Incubate EC (Fecal Coliform) tubes at $44.5 \pm 0.2^{\circ}$ C for 24 ± 2 h.
- Interpretation:
- Formation of gas in any amount, in Durham tube, of the brilliant green tubes, at any time, within 48 ± 3 h constitutes a positive confirmed result.
- Formation of gas in any amount, in Durham tube, of the EC tubes, at any time, within 24 ± 2 h constitutes a positive confirmed result.
- Calculate the MPN value of the number of positive brilliant green lactose bile tubes from MPN index (see A.3).

4. Completed Phase:

• Application:

- This phase is aiming to verify the presence of coliform bacteria and to provide quality control data for non-potable water sample analysis.
- Use the completed test on 10% of positive confirmed tubes on seasonal basis. Analysis of samples for thermotolerant (fecal) coliforms at elevated temperature (44.5±0.2°C) can be considered as a completed test.
- Alternatively, the completed test for positive total coliforms may be performed as follows:

• Media Preparation:

- Use LES Endo agar plates and Nutrient agar slants for completed phase.
- m-Endo LES agar:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense in 5 to 7 ml quantities into lower section of 50 mm Petri dishes, if dishes of any other size are used, adjust quantity to give an equivalent depth of 4 to 5 mm.
- If plates are previously prepared and stored in the refrigerator, allow them to warm at room temperature. The crystals that form on agar plates after refrigeration will disappear as the plates warm up. Use plates within two weeks after preparation.
- Nutrient agar:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml in screw-capped tubes before sterilization.
- Place tubes in an inclined position immediately after sterilization, so that the agar will solidify with a sloped surface.
- Tighten screw caps after cooling and store in a protected, cool storage area.
- If plates or tubes are made ahead of time and stored in the refrigerator, remove them and allow warming at room temperature before use.

• Procedure:

- Use a sterile 3-mm-diam loop or an inoculating needle slightly curved at the tip
- Tap and incline the fermentation tube to avoid picking up any membrane or scum on the needle
- Insert the end of the loop or needle into the liquid in the tube to a depth of approximately 0.5 cm
- Streak a plate for isolation with the curved section of the needle in contact with the agar to avoid a scratched or torn surface
- Change loop between second and third quadrants to improve colony isolation.
- Incubate plates (inverted) at $35 \pm 0.5^{\circ}$ C for 24 ± 2 h.
- Pick from each plate one or more typical (pink to dark red with a green metallic surface sheen) well-isolated coliform colonies. If no typical colonies are present, pick two or more colonies considered most likely to consist of organisms of the coliform group.

- Transfer growth from each isolate, by barely touching the surface of colony to minimize the danger of transferring mixed culture, to a single-strength lauryl tryptose broth fermentation tube and onto a nutrient agar slant.
- Incubate nutrient agar slants at 35 ± 0.5 °C for 24 ± 2 h.
- Incubate secondary broth tubes at 35 ± 0.5°C for 24 ± 2 h; if gas is not produced within 24 ± 2 h re-incubate and examine again at 48 ± 3 h.
- Examine Gram stained preparations from nutrient agar slant cultures corresponding to the secondary tubes that showed gas; include gram positive and gram-negative cultures as controls.

• Gram Stain technique:

- Prepare separate light emulsions of the test bacterial growth using drops of distilled water on the slide.
- Prepare separate light emulsions of the positive and negative control cultures on the same slide using drops of distilled water on the slide.
- Fix the slides by air drying.
- Stain for 1 min with ammonium oxalate-crystal violet solution, rinse slides in tap water and drain off excess.
- Apply Lugol's solution for 1 min, rinse stained slides in tap water.
- Decolorize for approximately 15 to 30 s with acetone alcohol by holding slide between the fingers and letting acetone alcohol flow across the stained smear until the solvent flows colorlessly from the slide. Do not over-decolorize to not remove crystal violet from gram positive bacteria.
- Counterstain with safranin for 15 s, rinse with tap water.
- Blot dry with absorbent paper or air dry, and examine microscopically.
- Gram-positive organisms are blue (violet); gram-negative organisms are red.
- Results are only acceptable when controls given proper reactions.



• Formation of gas in the secondary tube of lauryl tryptose broth within 48 ± 3 h and demonstration of gram-negative, nonspore-forming, rod-shaped bacteria from the agar culture constitute a positive result for completed test, demonstrating the presence of a member of the coliform group.

10. Calculations:

- The MPN values, for variety of positive and negative tubes combinations, are given in appendix A.3.
- The samples volumes indicated in indexes A.3 illustrates MPN values for combinations of positive results when five 10 ml, five 1 ml, and five 0.1 ml sample portion volumes are tested.
- If the sample portion volumes used are those found in the tables, report the value corresponding to the number of positive and negative results in the series as the MPN/100ml.
- When the series of decimal dilutions is different from that in the table, select the MPN value from index A.1 for the combination of positive results and calculate according to the following formula:

• MPN / 100 ml = (Table MPN/ 100 ml) \times 10/V

- Where:
 - V = volume of sample portion at the lowest selected dilution.
- When more than three dilutions are used in a decimal series of dilutions, select the three most appropriate dilutions and refer to index A.3.

• Selected examples(A though C) illustrating correct selection are listed below, more examples are available in Standard Methods:

Example	Volumes <i>ml</i>				Combination of	MPN index /	
Example	10	1.0	0.1	0.01	0.001	positive	100 ml
А	5	5	1	0	0	5-1-0	330
В	4	5	1	0	0	4-5-1	48
С	4	3	0	1	1	4-3-2	39

- Example A: First, remove the highest dilution (smallest sample volume) if it has all negative tubes and at least one remaining dilution has a negative tube. Next, remove the lowest dilution (largest sample volume) if it has all positive tube and at least one remaining dilution has a positive tube. According to these guidelines, the three dilutions in Example A are selected by removal of the highest (0.001-ml) and the lowest (10-ml) dilutions.
- Example B: If the lowest dilution does not have all positive tubes, and several of the highest dilutions have all negative tubes, then remove the highest negative dilutions.
- Example C: If no dilution has all positive tubes, select the lowest two dilutions, corresponding to 10 and 1 ml sample. For the third dilution, add the number of positive tubes in the remaining dilutions (0.1, 0.01, and 0.001 ml sample), to yield a final combination of 4-3-2
- If MPN value for combinations not appearing in the table, or for other combinations of tubes or dilutions, estimate the result as follows:
- First, select the lowest dilution that does not have all positive results.
- Second, select the highest dilution with at least one positive result.
- Select all dilutions between them.
- Use the selected dilutions in the following formula:

• MPN / 100 ml (approx.) = $100 \times P / (N \times T) \frac{1}{2}$

- Where:
 - P: number of positive results.
 - N: Volume of sample in all negative portions combined.
 - T: Total volume of sample in the selected dilutions.
- Examples:

- From (5/5, 10/10, 4/10, 1/10, 0/5) Use only (-, -, 4/10, 1/10.), 4/10 @ 0.1 ml sample/tube and 1/10 @ 0.01 ml sample/tube. calculations will be:
- MPN / 100 ml (approx.) = 100×5 / (0.69 ×1.1) $\frac{1}{2}$ = 500/0.87 = 570/ 100 ml
- From (5/5, 10/10, 10/10, 0/10, 0/5) Use only (-, -, 10/10, 0/10. -) 10/10 @ 0.1 ml sample/tube and 0/10 @ 0.01 ml sample/tube.
- calculations will be:
- MPN / 100 ml (approx.) = 100×10 / (0.1×1.1) $\frac{1}{2}$ = 1000/0.332 = 3000/100 ml.

11. Quality Control:

- Apply the presumptive-confirmed phase of the multiple-tube procedure to all samples, if samples were analyzed for fecal coliform using MTF technique there will be no need to apply completed phase, otherwise proceed completed test (Section 10.6) to not less than 10% of all coliform-positive samples on a seasonal basis.
- Plates with non-selective agar medium, labeled as "Environmental Control", will be opened along the samples analysis run time to ensure adequate environmental conditions.
- Procedural Blank: Inoculate, under the same conditions as sample, one additional tube for each dilution with sterile phosphate buffered dilution water after inoculation of series of 10 samples. Absence of growth indicates absence of cross contamination.

12. Reporting:

- Report coliform concentration as the Most Probable Number (MPN)/100 ml.
- If sample was delayed more than permissible time, write "Delayed" in report comment.

13. References:

- Eugene W. Rice, Rodger B. Baired, Andrew D. Eaton, Lenore S. Clesceri eds.2012 Standard Methods for the Examination of Water and Wastewater. Method 9221B#, 22nd edition (on-line edition). American Public Health Association, American Water Works Association, and Water Environment Federation
- World Health Organization. 2011, Guidelines for Drinking- Water Quality, Fourth Edition. Geneva
- Guidelines for Canadian Drinking water: Guideline Technical Document. 2006.
 Total Coliform. Health Canada. Ottawa. Ontario.

3. Detection and Enumeration of Total Coliform by Membrane Filter Technique

1. Scope:

- This is a detailed procedure for the detection and enumeration of Total Coliform (TC) by Membrane Filter (MF) technique in water samples in 24 hours or less on the basis of the production of aldehydes from fermentation of lactose.
- This procedure can be applied for different types of water, in microbiology laboratory it will be applied for the analysis of ground water, network distribution system (drinking water), and water from different treatment process.



2. Principle:

- Volume of water to be tested is filtered through 0.45 µm and the membrane is placed on Endo type medium. The selectivity of medium is due to sodium lauryl sulphate and sodium deoxycholate that acting as inhibitor of gram positive bacteria. Coliform ferment lactose, produce acetaldehydes that reacts with sodium sulfite / basic fuchsin compound to form red colony. The development of metallic sheen occurs when the organism produce aldehyde with the rapid fermentation of lactose. If the inoculum is too heavy, the sheen will be suppressed. Lactose non-fermenting bacteria form clear, colorless colonies.
- Total coliform bacteria (excluding E. coli) occur in both sewage and natural waters. Some of these bacteria are excreted in the faeces of humans and animals, but many coliforms are heterotrophic and able to multiply in water and soil environments.

• Total coliforms include organisms that can survive and grow in water distribution system, particularly in the presence of biofilms. They were, traditionally, regarded as belonging to the genera Escherichia, Citrobacter, Klebsiella and Enterobacter, but the group is more heterogeneous and includes a wider range of genera, such as Serratia and Hafnia. Hence, they are not useful as an index of faecal pathogens, but they can be used as an indicator of treatment effectiveness and to assess the cleanliness and integrity of distribution systems and the potential presence of biofilms. However, there are better indicators for these purposes.

Pre 1994 Acid and Gas from Lactose	Report 71, 1994 Acid from Lactose	Enzyme-based ß-Galactosidase
Escherichia	Escherichia	Escherichia
Klebsiella	Klebsiella	Klebsiella
Enterobacter	Enterobacter	Enterobacter
Citrobacter	Citrobacter	Citrobacter
	Yersinia	Yersinia
	Serratia	Serratia
	Hafnia	Hafnia
	Pantoea	Pantoea
	Kluyvera	Kluyvera
		Cedecea
		Ewingella
		Moellerella
		Leclercia
		Rahnella
		Yokenella

bold type = coliforms which can be present in the environment as well as in human faeces. **bold and underline** = coliforms which are considered to be primarily environmental.

Total Coliform species according to testing method

- The presence of Total Coliform after disinfection indicates inadequate treatment, and or stored water supplies can reveal regrowth and possible biofilm formation or contamination through ingress of foreign material, including soil or plants.
- Membrane Filter technique is highly reproducible technique which depends on sample filtration through a 47-mm, 0.45µm pore size cellulose membrane filters that retains the bacteria present in the sample. It can be used to test relatively large sample volumes and usually yields numerical results more than multiple-tube fermentation procedure.

3. Definitions:

• Total coliform bacteria in this method are those facultative anaerobic, gramnegative, non-spore-forming, rod-shaped bacteria that produce red colony with metallic (golden) sheen within 24 h incubation at 35°C on an Endo-type medium containing lactose. The sheen may cover the entire colony or may appear only in a central area or on the periphery.

4. Interference:

- Non-coliform bacteria may interfere with the recovery of coliforms when using a lactose-based medium. Data showed that the recovery of total coliforms using the MF technique decreased as the concentration of HPC bacteria increased. The greatest reduction occurred when the HPC densities exceeded 500 colony-forming units (CFU/ml). It should be noted that most water supplies maintaining a total chlorine residual of 0.2 mg/L have an HPC below 500 CFU/ml.
- Another data demonstrated that Pseudomonas aeruginosa (30 CFU/ml) and Aeromonas hydrophila (2 CFU/ml) caused significant reductions in sheen production by coliforms on m-Endo LES agar. Flavobacterium sp. and Bacillus sp., in contrast, were not inhibitory, even at concentrations above 1000 CFU/ml.
- Samples with high turbidity caused by algae, particulate, or other interfering material can clog the membrane filter, thereby preventing filtration, or cause spreading of bacterial colonies which could interfere with identification of target colonies.
- Some members of the total coliform group may produce dark red, mucoid, or nucleated colonies without a metallic sheen. When verified these are classified as atypical coliform colonies.
- Pink (non-mucoid), blue, white, or colorless colonies lacking sheen are considered non-coliforms by this technique.
- Also, do not use the MF technique to test wastewater containing high levels of toxic metals or toxic organic compounds (e.g. phenol) because such substances may be concentrated by the filter and inhibit coliform growth.

5. Equipment:

- Stainless steel forceps.
- Sterile plastic or glass petri dishes 50 mm.
- Sterile membrane filtration units.
- Vacuum pump.
- Buffered rinse water screw cap bottles.
- Measuring 100 ml Cylinder class A.
- Sterile membrane filters (47 mm diameter, 0.45 $\mu \pm 0.02 \mu m$ pore size, white, grid-marked).

- Sterile 47-mm diameter absorbent pads (used with broth).
- Electric Gas Burner.
- Incubator: $35^{\circ}C \pm 0.5^{\circ}C$.
- Colonies counter aid.
- Indelible ink marker for labeling plates.
- Sterile plastic loop, at least 3 mm diameter in suitable holders.

6. Chemicals and Reagents:

- M-Endo LES medium.
- Lauryl Tryptose broth.
- Brilliant Green Lactose Bile Broth.
- Phosphate buffered rinse water.
- Ethyl alcohol 95% in small wide-mouth vials.

7. Precautions:

- Follow the normal safety procedures required in a microbiology laboratory.
- Mouth-pipetting is prohibited.
- Decontaminate all the used plates and materials at the end of the analyses.
- Use sterile filtration units at the beginning of each filtration series as a minimum precaution to avoid accidental contamination. A filtration series is considered to be interrupted when an interval of 30 min or longer elapses between samples filtrations.
- After such interruption, treat any further samples filtration as a new filtration series and sterilize all membrane filter holders in use.
- Avoid damaging or dislodging the membrane filter when attaching funnel to base of the filter unit.
- Do not expose poured plates with culture media to direct light; refrigerate in the dark, in a suitable tight container.

8. Procedure:

1. Sample Handling:

- Analyze samples on day of receipt whenever possible. The time between sample collection and the placement of samples in the incubator must not exceed 30 hours.
- If arrival is too late for processing on same day, refrigerate overnight as long as holding time conditions still be met i.e. preservation temperature <10.0°C and time did not exceeded 30h.

2. Selection of sample size:

- The official volume for regulation purpose to be filtered is 100 ml.
- Use sample volumes that will yield counts between 20 and 80 total coliform colonies per membrane. Recommended volumes will be 100 and 50 ml.
- Analyze waters by filtering desired volume in same funnel, or by filtering replicate smaller sample volumes e.g. duplicate 50-ml or four replicates of 25-ml portions.
- For special monitoring purposes, such as troubleshooting water quality problems or identification of coliform breakthrough in low concentrations from treatment barriers; it can be test up to 1-L samples. If particulates prevent filtering a 1-L sample through a single filter, divide sample into four portions of 250 ml for analysis. Total the coliform counts on each membrane to report the number of coliforms per 100 ml.

3. Media and Reagents preparation:

• Prepare and handle media and reagents appropriately, carry out suitable performance test to ensure sterility and suitability of each prepared batch.

• Phosphate buffered rinse water:

- See annex No. A.1.
- m-Endo LES medium:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense in 5 to 7 ml quantities into lower section of 50 mm Petri dishes, if dishes of any other size are used, adjust quantity to give an equivalent depth of 4 to 5 mm.
- If liquid medium is used, place a pad in the culture dish and saturate with about 2.0 ml broth medium and carefully remove excess medium by decanting the plate.
- If plates are previously prepared and stored in the refrigerator, allow them to warm at room temperature. The crystals that form on agar plates after refrigeration will disappear as the plates warm up.
- Use plates within two weeks after preparation.
- Lauryl tryptose broth:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml medium, before sterilization, in fermentation tubes with an inverted Durham tube.
- Alternatively, and only in case of unavailability of Durham tubes, add 0.01 g/L bromcresol purple to medium to determine acid production.
- Close tubes with heat-resistant plastic caps.

- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes) after preparation.

• Brilliant green broth:

- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml medium, before sterilization, in fermentation tubes with an inverted Durham tube.
- Close tubes with heat-resistant plastic caps.
- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes) after preparation.

4. Sample Filtration:

- Label the bottom of plates with the sample identification number, analyst initials, and sample volume to be analyzed.
- Place a membrane filter (grid side up) using sterile forceps on the porous plate of the filter base.
- If you have difficulties in removing the separation papers from the filters due to static electricity, place a filter with the paper on top of the funnel base and turn on the vacuum, the separation paper will curl up, allowing easier removal.
- Attach the funnel to the base of the filter unit
- Shake the sample container vigorously.
- Measure the selected volume by using sterile measuring cylinder.
- Pour selected volume or dilution of the sample into the funnel.
- Turn on the vacuum pump and allow the vacuum to operate for a time enough to filtrate all the sample volume.
- Rinse the interior surface of funnel, with filter still in place, by filtering 20 to 30 ml portions of sterile dilution water. Alternatively, rinse funnel by a flow of sterile dilution water from a squeeze bottle.
- Turn off the pump then remove the funnel from the base of the filter unit.

- Hold the membrane filter at its edge with sterile forceps, gently lift and place the filter grid-side up on the medium plate.
- Slide the filter onto the medium, using a rolling action to avoid trapping air bubbles between the membrane filter and the underlying medium.
- Reseat the membrane if non-wetted areas occur due to air bubbles.

5. Incubation:

• Invert the plates, and incubate at 35 ± 0.5 °C for 22 to 24 h.



6. Counting:

- Determine colony count on membrane filters by using a suitable optical device with a cool white fluorescent light source to provide optimal viewing of sheen.
- Count plates that have 20 to 80 coliform colonies and not more than 200 colonies of all types on a membrane-filter surface.
- The typical coliform colony has a pink to dark-red color with a metallic surface sheen.
- Atypical coliform colonies can be dark red, mucoid, or nucleated without sheen.
- Pink, blue, white, or colorless colonies lacking sheen are considered non-coliforms.
- Count both typical and atypical coliform colonies.
- Refrigerate cultures (after incubation) with high densities of non-coliforms colonies for 0.5 to 1 h before counting to deter spread of confluence and aiding sheen discernment.

7. Verification:

• Preferably, verify all typical and atypical colony types,

but at a minimum, verify at least five typical and five atypical colonies from a given





membrane filter culture. If counts are > 80, verify by swabbing all colonies on membrane.

- Adjust counts on the basis of verification results.
- Verification Test:
 - Lactose Fermentation technique:
- Pick up the selected colony(ies) and inoculate simultaneously into both Lauryl tryptose broth tube(s) and Brilliant Green broth tube(s).
- Incubate all tubes at 35.0 ± 0.5 °C for 24-48 hr.
- Colony that revealed positive reaction in Lauryl tryptose broth tube (gas formed associated with growth turbidity, or color change to yellow in case of using Bromcresol purple instead of Durham tube) and in brilliant Green tube (gas formation) within 24-48 hr. is considered total coliform.

8. Calculations:

• Compute the count, using membrane filters with 20 to 80 coliform colonies and not more than 200 colonies of all types per membrane (see section 14.0), by the following equation:

Total Coliform / 100 mL = Coliforms colonies counted X 100 mL sample filtered

• For verified typical and atypical coliform counts, adjust the initial count based upon the positive verification percentage, sum both types, and report as "verified coliform count/100 ml." by using equation:

No of verified coloniesX Total no. of coloniesNo. of colonies subjected to verificationX Total no. of colonies

9. Quality Control:

- Filter each sample volume in duplicate manner, final value for that volume will be obtained by getting the arithmetic mean for both plates.
- Plate(s) with non-selective agar medium, labeled as "Environmental Control", will be opened along the samples analysis run time to ensure adequate environmental conditions.
- Filtration control (Procedural Blank): Filter 100 ml of sterile phosphate buffered rinse water, using separate sterile funnel, after filtration of series of 10 samples;

incubate under the same conditions as a sample. Absence of growth indicates the absence of cross- contamination.

10. Reporting:

- Report results as "Total Coliform / 100 ml of sample".
- If HPC result was >500 CFU/ml, report as "Estimated" in case of positive total coliform, or "False Negative" in case of negative result.
- If sample was delayed more than permissible time, write "Delayed" in report comment.
- If no coliform colonies are observed, report the coliform colonies counted as "<1 coliform/100 ml."
- If the total number of bacterial colonies, coliforms plus non-coliforms, exceeds 200 per membrane, report results as "Too Numerous To Count" (TNTC)
- If confluent growth occurs, covering either the entire filtration area of the membrane or a portion thereof, and colonies are not distinct enough for accurate counting, report results as "confluent growth with (or without) coliforms."
- Report confluent growth or TNTC without detectable coliform as 'Invalid''
- If sample was divided into two or more portions, Total the coliform counts on all filters and report the number of coliforms per 100 ml.
- Report the largest volume filter that has a coliform count falling in the ideal range; calculate final concentration value by multiplying the count by dilution factor (if present).
- If all volumes filters have a coliform count lower than the ideal range, disregard the rule and report the result from largest volume filter count; calculate final concentration value by multiplying the count by dilution factor (if present).
- If largest volume filter has a coliform count higher than the ideal range, choose next (second) dilution volume that falls inside ideal range and calculate final concentration value by multiplying the count by dilution factor.
- If next (second) dilution volume has coliform count lower than ideal range, and largest volume count is less than 100, use the count from the largest volume rather than the count from second volume; calculate final concentration value by multiplying the count by dilution factor.

• If next (second) dilution volume has coliform count higher than ideal range, report filter count and calculate final concentration value by multiplying the count by dilution factor.

1		Framn	عما	•
	••	слатир	103	•

Examples	100 ml	50 ml	Reported Value
01	54	18	54
02	10	5	10
03	89	52	104
04	90	18	90
05	101	18	36
06	TNTC	90	180
07	TNTC	TNTC	TNTC

الاداره العامه للمسار الوظيفي بالشركة القابضة لمياه الشرب والصرف الصحي

1. Detection and Enumeration of Thermotolerant (Fecal) Coliform by Membrane Filter Technique

1. Scope:

- This method describes a detailed procedure for detection and enumeration of Thermotolerant coliform (Fecal Coliform) by Membrane Filter (MF) technique in water samples in 24 hours on the basis of fermentation of lactose at elevated temperature.
- This procedure can be applied for different types of water: (ground water, network distribution system (drinking water), and water from different treatment process).



2. Principle:

- Volume of water to be tested is filtered through 0.45 μm and the membrane is placed on m-FC medium.
- The selectivity of medium is due to Bile salt No.3 that acting as inhibitor of gram positive bacteria in addition to elevated incubation temperature (44.5± 0.2°C). Aniline blue is pH indicator turning into blue color in acidic medium.
- Thermotolerant coliform (Fecal Coliform) ferment lactose, produce acid that change pH of the medium into acidic, which results in coloring colonies with blue.
- Total coliform bacteria that are able to ferment lactose at 44.5 °C are known as thermo-tolerant coliforms. Thermo-tolerant coliforms were traditionally called fecal coliforms, but they also have been documented in organically rich waters or tropical climates in the absence of recent fecal contamination. So, testing for E. coli, a specific indicator of fecal contamination, is recommended.

- In most waters, the predominant genus is Escherichia, but some types of Citrobacter, Klebsiella and Enterobacter are also thermotolerant. They are usually found in sewage and water recently subjected to fecal pollution.
- Populations of thermo-tolerant coliforms are composed predominantly of E. coli; as a result, this group is regarded as a less reliable but acceptable index of faecal pollution. The presence of thermo-tolerant coliform (fecal coliform) provides evidence of recent fecal pollution.
- Membrane Filter technique is highly reproducible technique which depends on sample filtration through a 47-mm, 0.45 µm pore size cellulose membrane filters that retains the bacteria present in the sample. It can be used to test relatively large sample volumes and usually yields numerical results more than multiple-tube fermentation procedure.

3. Definitions:

• Thermotolerant coliform (Fecal Coliform) bacteria in this method are those facultative anaerobic, gram-negative, non-spore-forming, rod-shaped bacteria that produce blue colony within 24 h incubation at $44.5 \pm 0.2^{\circ}$ C on m FC medium.

4. Interference:

- Samples of storm water collected during the first runoff (initial flushing) after a long dry period may have a background growth that will interfere with the recovery of fecal coliform. To eliminate such interference, add 1 % rosolic acid salt reagent to the prepared media.
- Samples with high turbidity caused by algae, particulate, or other interfering material can clog the membrane filter, thereby preventing filtration, or cause spreading of bacterial colonies which could interfere with identification of target colonies.
- Non-fecal coliform colonies, a gray to cream-colored, may be observed on m FC media.

5. Equipment:

- Stainless steel forceps.
- Sterile plastic petri dishes 50 mm.
- Sterile membrane filtration units.
- Vacuum pump.

- Buffered rinse water screwed cap bottles.
- Measuring 100 ml Cylinder class A.
- Sterile membrane filters (47 mm diameter, 0.45μ ± 0.02 μm pore size, white, grid-marked).
- Sterile 47-mm diameter absorbent pads (used with broth).
- Electric Gas Burner.
- Incubator: $44.5^{\circ}C \pm 0.2^{\circ}C$.
- Incubator: $35^{\circ}C \pm 0.5^{\circ}C$.
- Colonies counter aid.
- Indelible ink marker for labeling plates.
- Sterile plastic loop, at least 3 mm diameter in suitable holders.

6. Chemicals and Reagents:

- m-FC medium.
- Lauryl Tryptose broth.
- EC broth.
- Phosphate buffered rinse water.
- Ethyl alcohol 95% in small wide-mouth vials.

7. Precautions:

- Follow the normal safety procedures required in a microbiology laboratory.
- Mouth-pipetting is prohibited.
- Decontaminate all the used plates and materials at the end of the analyses.
- Use sterile filtration units at the beginning of each filtration series as a minimum precaution to avoid accidental contamination. A filtration series is considered to be interrupted when an interval of 30 min or longer elapses between samples filtrations.
- After such interruption, treat any further samples filtration as a new filtration series and sterilize all membrane filter holders in use.
- Avoid damaging or dislodging the membrane filter when attaching funnel to base of the filter unit.
- Do not expose poured plates with culture media to direct light; refrigerate in the dark, in a suitable tight container.

8. Procedure:

1. Sample Handling:

- Analyze samples on day of receipt whenever possible. The time between sample collection and the placement of samples in the incubator must not exceed 30 hours.
- If arrival is too late for processing on same day, refrigerate overnight as long as holding time conditions still be met i.e. preservation temperature <8.0°C and time did not exceeded 30 hr.

2. Selection of sample size:

- The official volume for regulation purpose to be filtered is 100 ml.
- Use sample volumes that will yield counts between 20 and 60 fecal coliform colonies per membrane. Recommended volumes will be 100 and 50 ml.
- Analyze waters by filtering desired volume in same funnel, or by filtering replicate smaller sample volumes e.g. duplicate 50-ml or four replicates of 25-ml portions.
- For special monitoring purposes, such as troubleshooting water quality problems or identification of fecal coliform breakthrough in low concentrations from treatment barriers; it can be test up to 1-L samples. If particulates prevent filtering a 1-L sample through a single filter, divide sample into four portions of 250 ml for analysis. Total the fecal coliform counts on each membrane to report the number of fecal coliforms per 100 ml.

3. Media and solutions preparation:

• Prepare and handle media and reagents appropriately, and carry out suitable performance test to ensure homogeneity, sterility, and suitability for each prepared batch.

• Phosphate buffered rinse water:

- See annex No. A.1.
- m-FC medium:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense in 5 to 7 ml quantities into lower section of 50 mm Petri dishes, if dishes of any other size are used, adjust quantity to give an equivalent depth of 4 to 5 mm.
- If liquid medium is used, place a pad in the culture dish and saturate with about 2.0 ml broth medium and carefully remove excess medium by decanting the plate.

- If plates are previously prepared and stored in the refrigerator, allow them to warm at room temperature. The crystals that form on agar plates after refrigeration will disappear as the plates warm up.
- Use plates within two weeks after preparation.
- Lauryl tryptose broth:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml medium, before sterilization, in fermentation tubes with an inverted Durham tube.
- Alternatively, and only in case of unavailability of Durham tubes, add 0.01 g/L bromcresol purple to medium to determine acid production.
- Close tubes with heat-resistant plastic caps.
- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes).
- EC broth:
- Prepare and sterilize medium according to manufacturer instructions.
- Dispense 10 ml medium, before sterilization, in fermentation tubes with an inverted Durham tube.
- Close tubes with heat-resistant plastic caps.
- Ensure, after sterilization, that media cover inverted tube with at least one-half to two-thirds, and completely filled with media i.e. no air space inside inverted Durham tubes.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes).

4. Sample Filtration:

- Label the bottom of plates with the sample identification, analyst initials, and sample volume to be analyzed.
- Place a membrane filter (grid side up) using sterile forceps on the porous plate of the filter base.

- If you have difficulties in removing the separation papers from the filters due to static electricity, place a filter with the paper on top of the funnel base and turn on the vacuum, the separation paper will curl up, allowing easier removal.
- Attach the funnel to the base of the filter unit
- Shake the sample container vigorously.
- Measure the selected volume by using sterile measuring cylinder.
- Pour selected volume or dilution of the sample into the funnel.
- Turn on the vacuum pump and allow the vacuum to operate for a time enough to filtrate all the sample volume.
- Rinse the interior surface of funnel, with filter still in place, by filtering three 20 to 30 ml portions of sterile dilution water. Alternatively, rinse funnel by a flow of sterile dilution water from a squeeze bottle.
- Turn off the pump then remove the funnel from the base of the filter unit.
- Hold the membrane filter at its edge with a sterile forceps, gently lift and place the filter grid-side up on the medium plate.
- Slide the filter onto the medium, using a rolling action to avoid trapping air bubbles between the membrane filter and the underlying medium.
- Reseat the membrane if non-wetted areas occur due to air bubbles.

5. Incubation:

• Invert the plates, and incubate at 44.5 ± 0.2 °C for 24 ± 2 h.



6. Counting:

- To determine colony count on membrane filters, use a suitable optical device with a cool white fluorescent light source to provide optimal viewing of sheen.
- Count colonies with various shades of blue as typical fecal coliform, grey to green colonies as atypical fecal coliform. Non-fecal coliform colonies are gray to cream-colored.

- Count both typical and a typical fecal coliform colony.
- 7. Verification:





- Verify at a frequency established by the laboratory. Verify typical blue colonies and any atypical grey to
- green colonies. (preferably verify all colonies per membrane)
- Adjust counts on the basis of verification results.
- Verification Test:
- Lactose Fermentation technique:
- Pick up the selected colony(ies) and inoculate simultaneously into both Lauryl tryptose broth tube(s) and EC tube(s).
- Incubate Lauryl tryptose broth tube(s) at 35.0 ± 0.5 °C for 24-48 hr.
- Incubate EC tube(s) at 44.5 ± 0.2 °C for 24 hr.
- Colony that revealed positive reaction in Lauryl tryptose broth tube (gas formed associated with growth turbidity, or color change to yellow in case of using Bromcresol purple instead of Durham tube) within 48 hr and positive reaction in EC tube (gas formation) within 24 hr is considered fecal coliform.

8. Calculations:

• Compute the count, using membrane filters with 20 to 60 fecal coliform colonies by the following equation:

Fecal coliform / 100 mL= Colonies counted × 100 ML sample filtered

• For verified typical and atypical fecal coliform counts, adjust the initial count based upon the positive verification percentage, sum both types, and report as "verified fecal coliform count/100 ml."

No. of verified colonies

X Total no. of colonies

No. of colonies subjected to verification

الاداره العامه للمسار الوظيفي بالشركة القابضة لمياه الشرب والصرف الصحي

9. Quality Control:

- Filter each sample volume in duplicate manner, final value for that volume will be obtained by getting the arithmetic mean for both plates.
- Plate(s) with non-selective agar medium, labeled as "Environmental Control", will be opened along the samples analysis run time to ensure adequate environmental conditions.
- Filtration control (Procedural Blank): Filter 100 ml of sterile phosphate buffered rinse water, using separate sterile funnel, after filtration of series of 10 samples; incubate under the same conditions as a sample. Absence of growth indicates absence of cross-contamination.

10. Reporting:

- Report results as "Fecal Coliform / 100 ml of sample"
- If sample was delayed more than permissible time, write "Delayed" in report comment.
- If no thermotolerant coliform (fecal coliform) colonies are observed, report the fecal coliform colonies as ''<1 fecal coliform/100 ml.''
- If the total number of thermotolerant coliform (fecal coliform) colonies exceeds 60 per membrane, report results as "Too Numerous To Count with (or without) fecal coliforms" (TNTC).
- If confluent growth occurs, covering either the entire filtration area of the membrane or a portion thereof, and colonies are not distinct enough for accurate counting, report results as "confluent growth with (or without) fecal coliforms."
- Report confluent growth or TNTC without detectable fecal coliform as 'Invalid''.
- If sample was divided into two or more portions, Total the fecal coliform counts on all filters and report the number of fecal coliforms per 100 ml.
- Report the largest volume filter that has a fecal coliform count falling in the ideal range; calculate final concentration value by multiplying the count by dilution factor (if present).
- If all volumes filters have a fecal coliform count lower than the ideal range, disregard the rule and report the result from largest volume filter count. Calculate final concentration value by multiplying the count by dilution factor (if present).

- If largest volume filter has a fecal coliform count higher than the ideal range, choose next (second) dilution volume that falls inside ideal range and calculate final concentration value by multiplying the count by dilution factor.
- If next (second) dilution volume has fecal coliform count lower than ideal range, and largest volume count is less than 100, use the count from the largest volume rather than the count from second volume; calculate final concentration value by multiplying the count by dilution factor.

Examples	100 ml	50 ml	Reported Value
01	54	18	54
02	10	5	10
03	89	52	104
04	90	18	90 or TNTC
05	101	18	36
06	TNTC	90	180
07	TNTC	TNTC	TNTC

2. Detection and Enumeration of Fecal Enterococcus /Streptococcus Groups by Membrane Filter Technique

1. Scope:

- This is a detailed procedure for detection and enumeration of Fecal Streptococcus group (FS) and Enterococcus group by Membrane Filter (MF) technique in water samples in 48 hours or less on the basis of reduction of Tri-phenyl Tetra-zolium Chloride (TTC).
- This procedure can be applied for different types of water: (ground water, network distribution system (drinking water), and water from different treatment process).

2. **Principle:**

 Volume of water to be tested is filtered through 0.45 µm and the membrane is placed on m-Enterococcus agar medium. The selectivity of medium is due to sodium azide that acts as inhibitor of gram negative bacteria. Bacteria reduce Triphenyl Tetrazolium Chloride (TTC) to the insoluble formazan inside the bacterial cell, resulting in the production of red colonies.



TCC is the dye used as an indicator of bacterial growth.

Fecal streptococcus colonies

- The fecal streptococcus group consists of a number of species of the genus Streptococcus, such as *S. bovis* and *S. equinus*. The normal habitat of fecal streptococci is the gastrointestinal tract of warm-blooded animals and humans.
- Enterococcus group includes *E. faecalis, E. faecium, E. gallinarum, and E. avium.* Enterococci group can be used as an index of fecal pollution, important advantages of this group are that they tend to survive longer in water environment than fecal coliforms, more resistant to drying and are more resistant to chlorination.

 Membrane Filter technique is highly reproducible technique which depends on sample filtration through a 47-mm, 0.45 or 0.22-µm pore size cellulose membrane filters that retains the bacteria present in the sample. It can be used to test relatively large sample volumes and usually yields numerical results.

3. Definitions:

Fecal Enterococcus bacteria in this method are those gram-positive, cocci-shaped, facultative anaerobic, catalase-negative, able to grow on bile esculin agar and in 6.5% NaCl broth at 35°C and either in BHI agar at 10 ± 0.5°C or BHI broth at 45 ± 0.5°C, producing red colony on m-Enterococcus agar.

4. Interference:

• Samples with high turbidity caused by algae, particulate, or other interfering material can clog the membrane filter, thereby preventing filtration, or cause spreading of bacterial colonies which could interfere with identification of target colonies.

5. Equipment:

- Stainless steel forceps.
- Presterilized plastic petri dishes 50 mm.
- Sterile membrane filtration units.
- Vacuum pump.
- Buffered rinse water screw cap bottles.
- Measuring 100 ml Cylinder class A.
- Sterile membrane filters (47 mm diameter, 0.45 $\mu \pm 0.02 \mu m$ pore size, white, grid-marked).
- Sterile 47-mm diameter absorbent pads (used with broth).
- Electric Gas Burner.
- Incubator: $35^{\circ}C \pm 0.5^{\circ}C$.
- Incubator: $45^{\circ}C \pm 0.5^{\circ}C$.
- Colonies counter aid.
- Indelible ink marker for labeling plates.
- Sterile plastic loop, at least 3 mm diameter in suitable holders.

6. Chemicals and Reagents:

- M-Enterococcus agar.
- Brain-heart infusion broth.

- Brain-heart infusion broth with 6.5 % NaCl.
- Brain-heart infusion agar.
- Bile esculin agar.
- Phosphate buffered rinse water.
- Hydrogen Peroxide 3.0 %
- Sodium Chloride (NaCl)
- Gram Stain reagents set.

7. Procedure:

1. Sample Handling:

- Analyze samples on a day of receipt whenever possible. The time between sample collection and the placement of samples in the incubator must not exceed 30 hours.
- If arrival is too late for processing on same day, refrigerate overnight as long as holding time conditions still be met i.e. preservation temperature <10.0°C and time did not exceeded 30h.

2. Selection of sample size:

- The official volume for regulation purpose to be filtered is 100 ml.
- Use sample volumes that will yield counts between 20 and 60 fecal streptococci colonies per membrane. Recommended volumes will be 100 and 50 ml.
- Analyze waters by filtering desired volume in same funnel, or by filtering replicate smaller sample volumes e.g. duplicate 50-ml or four replicates of 25-ml portions.
- For special monitoring purposes, such as troubleshooting water quality problems or identification of fecal streptococcus or enterococcus breakthrough in low concentrations from treatment barriers; it can be test up to 1-L samples. If particulates prevent filtering a 1-L sample through a single filter, divide sample into four portions of 250 ml for analysis. Total the counts on each membrane to report the number of fecal streptococcus or enterococcus per 100 ml.

3. Media and solutions preparation:

- Prepare and handle media and reagents appropriately, carry out suitable performance test to ensure sterility and suitability of each prepared batch.
- Phosphate buffered rinse water:
- See annex No. A.1.
- m-Enterococcus agar medium:
- Prepare and sterilize medium according to manufacturer instructions.

- Dispense in 5 to 7 ml quantities into lower section of 50 mm Petri dishes, if dishes of any other size are used, adjust quantity to give an equivalent depth of 4 to 5 mm.
- If plates are previously prepared and stored in the refrigerator, allow them to warm at room temperature. The crystals that form on agar plates after refrigeration will disappear as the plates warm up.
- Use plates within two weeks after preparation.
- WARNING The selective m-Enterococcus agar medium contain sodium azide. As this substance is highly toxic and mutagenic, precautions shall be taken to avoid contact with it, especially by the inhalation of fine dust during the preparation of commercially available dehydrated complete media. Azide-containing media should not be mixed with strong inorganic acids, as toxic hydrogen azide (HN3) may be produced. Solutions containing azide can also form explosive compounds when in contact with metal pipework, for example from sinks.

• Brain-heart infusion broth:

- Prepare medium according to manufacturer instructions.
- Dispense 10 ml medium in tubes, before sterilization.
- Close tubes with metal or heat-resistant plastic caps.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes).

• Brain-heart infusion broth with 6.5 % NaCl:

- Prepare medium according to manufacturer instructions.
- Add 60.0 g/l of NaCl to the ingredient of brain heart infusion broth.
- Close tubes with metal or heat-resistant plastic caps.
- Use tubes within two weeks (if loose plastic cap used) and within three months (for tight screw cap tubes) after preparation.

• Brain-heart infusion agar:

- Prepare medium according to manufacturer instructions.
- Add 15.0 g agar to the ingredient of brain heart infusion broth.
- Dispense 8 to 10 ml in screw-capped tubes, before sterilization.
- Place tubes in an inclined position, immediately after sterilization, so that the agar will solidify with a sloped surface.
- Tighten screw caps after cooling and store in a protected, cool storage area.
- Use tubes within three months after preparation.

• Bile esculin agar:

- Prepare medium according to manufacturer instructions
- Dispense 8 to 10 ml in screw-capped tubes, before sterilization.
- Place tubes in an inclined position, immediately after sterilization, so that the agar will solidify with a sloped surface.
- Tighten screw caps after cooling and store in a protected, cool storage area.
- Use tubes within three months after preparation.

4. Sample Filtration:

- Label the bottom of plates with the sample identification, analyst initials, and sample volume to be analyzed.
- Place a membrane filter (grid side up) using sterile forceps on the porous plate of the filter base.
- If you have difficulties in removing the separation papers from the filters due to static electricity, place a filter with the paper on top of the funnel base and turn on the vacuum, the separation paper will curl up, allowing easier removal.
- Attach the funnel to the base of the filter unit
- Shake the sample container vigorously.
- Measure the selected volume by using sterile measuring cylinder.
- Pour selected volume or dilution of the sample into the funnel.
- Turn on the vacuum pump and allow the vacuum to operate for a time enough to filtrate all the sample volume.
- Rinse the interior surface of funnel, with filter still in place, by filtering three 20 to 30 ml portions of sterile dilution water. Alternatively, rinse funnel by a flow of sterile dilution water from a squeeze bottle.
- Turn off the pump then remove the funnel from the base of the filter unit.
- Hold the membrane filter at its edge with sterile forceps, gently lift and place the filter grid-side up on the medium plate.
- Slide the filter onto the medium, using a rolling action to avoid trapping air bubbles between the membrane filter and the underlying medium.
- Reseat the membrane if non-wetted areas occur due to air bubbles.

5. Incubation:

• Let the plates stand for 30 min, invert the plates, and incubate at 35 ± 0.5 °C for 48 h

6. Counting:

- Determine colony count on membrane filters by using a suitable optical device with a cool white fluorescent light source to provide optimal viewing of sheen.
- Count plates that have 20 to 60.
- The typical colony has red color.
- Count all light and dark red colonies as presumptive enterococci.



Fecal streptococcus colonies

7. Verification:

- Because this method suffer periodically from false positives, verification test are an important quality-control step.
- Include a routine verification procedure, preferably similar to total coliform procedure i.e. at least five colonies from both dark and red colonies.
- Adjust counts on the basis of verification results.
- Verification Test: (See annex no. 3)
- Biochemical Characteristics:
- Pick selected typical colony from a membrane and streak for isolation onto the surface of a brain-heart infusion agar.
- Incubate at $35 \pm 0.5^{\circ}$ C for 24 to 48 h.
- Transfer a loop full of growth from a well-isolated colony on brain-heart infusion agar into a brain-heart infusion broth tube and to each of two clean glass slides.
- Incubate the brain-heart infusion broth at 35 ± 0.5 °C for 24 h.
- Add a few drops of freshly prepared 3% hydrogen peroxide to the smear on a slide. The appearance of bubbles constitutes a positive catalase test and indicates that the colony is not a member of the fecal streptococcus group. If the catalase test is negative, i.e., no bubbles, make a Gram stain of the second slide.
- Transfer a loop full of growth from the brain-heart infusion broth to each of the following media:

- bile esculin agar (incubate at 35 ± 0.5 °C for 48 h);
- brain-heart infusion agar (incubate at 10 ± 0.5°C for 48 h) or brain-heart infusion broth (incubate at 45 ± 0.5°C for 48 h);
- brain-heart infusion broth with 6.5% NaCl (incubate at 35 ± 0.5 °C for 48 h).
- Growth of catalase-negative, gram-positive cocci on bile esculin agar and in 6.5% NaCl broth at 35°C and either in BHI agar at 10 ± 0.5 °C or BHI broth at 45 ± 0.5 °C confirms that the colony belongs to the Enterococcus genus.

8. Calculations:

• Compute the count, using membrane filters with 20 to 60 colonies, by the following equation:

Fecal Streptococci/100 mL= ML sample filtered

• For verified fecal streptococci and enterococci counts, adjust the initial count based upon the positive verification percentage, sum both types, and report as "verified fecal streptococcus /100 ml." by using equation:

Number of verified colonies × 100Total number of colonies subjected to verification

8. Quality Control:

- Filter each sample volume in duplicate manner, final value for that volume will be obtained by getting the arithmetic mean for both plates.
- Plate(s) with non-selective agar medium, labeled as "Environmental Control", will be opened along the samples analysis run time to ensure adequate environmental conditions.
- Filtration control (Procedural Blank): Filter 100 ml of sterile phosphate buffered rinse water, using separate sterile funnel, after filtration of series of 10 samples; incubate under the same conditions as a sample. Absence of growth indicates absence of cross-contamination.

9. Reporting:

- Report the result as "Fecal enterococci / 100 ml"
- If sample was delayed more than permissible time, write "Delayed" in report comment.
- If no colonies are observed, report the fecal enterococci colonies counted as "<1 fecal enterococci /100 ml."
- If no filter has count falling in the ideal range, total the counts (disregarding the lower limit of 20 cited above) and use the formula given above to obtain fecal enterococci density.
- If confluent growth occurs, covering either the entire filtration area of the membrane or a portion thereof, and colonies are not distinct enough for accurate counting, report results as "confluent growth with (or without) fecal enterococci."
- If the total number of bacterial colonies exceeds 200 per membrane, report results as "Too Numerous To Count" (TNTC)
- If sample was divided into two or more portions, Total the counts on all filters and report the number of fecal enterococci per 100 ml.
- Report the largest volume filter that has count falling in the ideal range; calculate final concentration value by multiplying the count by dilution factor (if present).
- If all volumes filter has count lower than the ideal range, disregard the rule and report the result from largest volume filter count. Calculate final concentration value by multiplying the count by dilution factor (if present).
- If largest volume filter has count higher than the ideal range, choose next (second) dilution volume that falls inside ideal range and calculate final concentration value by multiplying the count by dilution factor.
- If next (second) dilution volume has count lower than ideal range, and largest volume count is less than 100, use the count from the largest volume rather than the count from second volume; calculate final concentration value by multiplying the count by dilution factor.
- If next (second) dilution volume has count higher than ideal range, report filter count and calculate final concentration value by multiplying the count by dilution factor.
- Examples:

Examples	100 ml	50 ml	Reported Value
01	54	18	54
02	10	5	10
03	89	52	104

04	90	18	90
05	101	18	36
06	TNTC	90	180
07	TNTC	TNTC	TNTC

10. References:

- A.D., Clesceri, L.S., Rice E. W. eds. 2005. Standard Methods for the Examination of Water and Wastewater, Method 9230 C# 21st Edition. American Public Health Association, American Water Works Association.
- U.S. Environmental Protection Agency.2005: Manual of The Certification of laboratories Analyzing Drinking Water, Criteria and Procedures Quality Assurance. EPA 815-R-05-054. Office of Ground Water and Drinking Water. Cincinnati, Ohio.
- U.S Environmental Protection Agency, Office of Research and Development, National Exposure research laboratory. 2003. Standard Operating Procedure for Choosing the best membrane filter count for the calculation of final concentration of Microorganisms per 100 ml. Number MERP-041.01

Annexes

A.1 Preparation of Phosphate Buffered Dilution Water

- Stock phosphate buffer solution: Dissolve 34.0 g KH₂PO₄ in a 500 ml reagentgrade distilled water, adjust the pH of the solution to 7.2 with 1 N NaOH and bring volume to 1000 ml with reagent-grade distilled water. Sterilize by filtration or autoclave for 15 minutes at 121°C, store in the refrigerator until using (If evidence of mold or other contamination appears in the stock, the solution should be discarded, and a fresh solution should be prepared).
- **MgCl₂ solution**: Dissolve 38 g anhydrous MgCl₂ (or 81.1 g MgCl₂.6H₂O) in one liter of reagent-grade distilled water. Sterilize by filtration or autoclave for 15 minutes at 121°C, store in the refrigerator until using (If evidence of mold or other contamination appears in the stock, the solution should be discarded, and a fresh solution should be prepared).
- Working solution: Add 1.25 ml phosphate buffer stock and 5 ml MgCl₂ stock for each liter of reagent-grade distilled water prepared, mix well(Final pH 7.0 ± 0.2), and dispense in appropriate amounts for dilutions in screw-cap dilution bottles or culture tubes, and/or into larger containers for use as rinse water. Autoclave at 121°C for 15 minutes.

A.2 Index: MPN Index and 95% confidence limits for concentrations of Positive and Negative results when Five 20-ml or Ten 10-ml portions are used.

No. of Tubes Giving Positive Reaction Out	MPN Index/	95% Co Limits	95% Confidence Limits (Exact)	
of 5 (20 mL Each)	100 mL	Lower	Upper	
0	<1.1	_	3.5	
1	1.1	0.051	5.4	
2	2.6	0.40	8.4	
3	4.6	1.0	13	
4	8.0	2.1	23	
5	> 8.0	3.4	-	

TABLE 9221:II. MPN INDEX AND 95% CONFIDENCE LIMITS FOR ALL COMBINATIONS OF POSITIVE AND NEGATIVE RESULTS WHEN FIVE 20-ML PORTIONS ARE USED

TABLE 9221:III. MPN INDEX AND 95% CONFIDENCE LIMITS FOR ALL COMBINATIONS OF POSITIVE AND NEGATIVE RESULTS WHEN TEN 10-ML PORTIONS ARE USED

No. of Tubes Giving Positive Reaction Out	MPN Index/	95% Co Limits	95% Confidence Limits (Exact)		
of 10 (10 mL Each)	100 mL	Lower	Upper		
0	<1.1	_	3.4		
1	1.1	0.051	5.9		
2	2.2	0.37	8.2		
3	3.6	0.91	9.7		
4	5.1	1.6	13		
5	6.9	2.5	15		
6	9.2	3.3	19		
7	12	4.8	24		
8	16	5.8	34		
9	23	8.1	53		
10	>23	13	-		

A.3 MPN index and 95% confidence limits for various combinations of positive results when five tubes are used per dilution (10 ml, 1.0 ml, 0.1 ml).

Combination of	Confidence Limits					Confidence Limits	
Positives	MPN Index/100 mL	Low	High	Positives	MPN Index/100 mL	Low	High
0-0-0	<1.8	_	6.8	4-0-3	25	9.8	70
0-0-1	1.8	0.090	6.8	4-1-0	17	6.0	40
0-1-0	1.8	0.090	6.9	4-1-1	21	6.8	42
0-1-1	3.6	0.70	10	4-1-2	26	9.8	70
0-2-0	3.7	0.70	10	4-1-3	31	10	70
0-2-1	5.5	1.8	15	4-2-0	22	6.8	50
0-3-0	5.6	1.8	15	4-2-1	26	9.8	70
1-0-0	2.0	0.10	10	4-2-2	32	10	70
1-0-1	4.0	0.70	10	4-2-3	38	14	100
1-0-2	6.0	1.8	15	4-3-0	27	9.9	70
1-1-0	4.0	0.71	12	4-3-1	33	10	70
1-1-1	6.1	1.8	15	4-3-2	39	14	100
1-1-2	8.1	3.4	22	4-4-0	34	14	100
1-2-0	6.1	1.8	15	4-4-1	40	14	100
1-2-1	8.2	3.4	22	4-4-2	47	15	120
1-3-0	8.3	3.4	22	4-5-0	41	14	100
1-3-1	10	3.5	22	4-5-1	48	15	120
1-4-0	11	3.5	22	5-0-0	23	6.8	70
2-0-0	4.5	0.79	15	5-0-1	31	10	70
2-0-1	6.8	1.8	15	5-0-2	43	14	100
2-0-2	9.1	3.4	22	5-0-3	58	22	150
2-1-0	6.8	1.8	17	5-1-0	33	10	100
2-1-1	9.2	3.4	22	5-1-1	46	14	120
2-1-2	12	4.1	26	5-1-2	63	22	150
2-2-0	9.3	3.4	22	5-1-3	84	34	220
2-2-1	12	4.1	26	5-2-0	49	15	150
2-2-2	14	5.9	36	5-2-1	70	22	170
2-3-0	12	4.1	26	5-2-2	94	34	230
2-3-1	14	5.9	36	5-2-3	120	36	250
2-4-0	15	5.9	36	5-2-4	150	58	400
3-0-0	7.8	2.1	22	5-3-0	79	22	220
3-0-1	11	3.5	23	5-3-1	110	34	250
3-0-2	13	5.6	35	5-3-2	140	52	400
3-1-0	11	3.5	26	5-3-3	170	70	400
3-1-1	14	5.6	36	5-3-4	210	70	400
3-1-2	17	6.0	36	5-4-0	130	36	400
3-2-0	14	5.7	36	5-4-1	170	58	400
3-2-1	17	6.8	40	5-4-2	220	70	440
3-2-2	20	6.8	40	5-4-3	280	100	710
3-3-0	17	6.8	40	5-4-4	350	100	710
3-3-1	21	6.8	40	5-4-5	430	150	1100
3-3-2	24	9.8	70	5-5-0	240	70	710
3-4-0	21	6.8	40	5-5-1	350	100	1100
3-4-1	24	9.8	70	5-5-2	540	150	1700
3-5-0	25	9.8	70	5-5-3	920	220	2600
4-0-0	13	4.1	35	5-5-4	1600	400	4600
4-0-1	17	5.9	36	5-5-5	>1600	700	_
4-0-2	21	6.8	40				

TABLE 9221:IV. MPN INDEX AND 95% CONFIDENCE LIMITS FOR VARIOUS COMBINATIONS OF POSITIVE RESULTS WHEN FIVE TUBES ARE USED FER DILUTION (10 Min, 1.0 Min, 0.1 ML)*

* Results to two significant figures.

A.4. Non-Conformance and Corrective Action:

- If the media performance and sterility assessment are performed concurrently with the use of the media in testing and the media quality is deemed unsatisfactory, the specific use of the media in the test(s) should be investigated to determine the potential impact on the test results and the necessity for repeating the test(s).
- Testing must be repeated when unsatisfactory media and/or reagents are used for critical aspects of testing.
- If a procedural blank test gives unsatisfactory results, the samples results should be investigated to determine the potential impact of contamination on the test result, if test results deemed unsatisfactory the test must be repeated.
- If, one or more of, environmental control plate(s) gives unsatisfactory results, the samples results should be investigated to determine the potential impact on the test result, if test results deemed unsatisfactory the test must be repeated.
- For all non-conformance findings, trace all media, reagents, equipment operations, and tools used in analysis to find out root cause of non-conformance.
- Take appropriate corrective action, if cause was temporary e.g. if unsterile batch of media used accidently in the analysis; discard all batch of prepared media, or preventive actions if cause was permanent e.g. re-maintenance of autoclave if actual temperature was less than adjusted value, to eliminate reoccurrence of non-conformance in the future.
- For Invalid results, request a new sample from the same location, explaining the reason for re-sampling need.

التحاليل البكتيريولوجية

A.3. Verification test for enterococcus group:



قام بإعداد الإصدار الثانى من هذا البرنامج:

كيميائى/ عاصم عبد الرحمن المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ أحمد كمال عبدالهادى المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ وائل عبدالرحيم أبوالمجد المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ كريم فاروق إسماعيل المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ محمود جمعة حسين المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ الحسن عبدالرافع عبدالصادق المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ الحسن عبدالرافع عبدالصادق المعمل المرجعى لمياه الشرب الشركة القابضة المعمل المرجعى لمياه الشرب الشركة القابضة المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ الحسن عبدالرافع عبدالصادق المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ الحسن عبدالرافع عبدالصادق المعمل المرجعى لمياه الشرب الشركة القابضة كيميائى/ الحسن عبدالرافع عبدالصادق المعمل المرجعى لمياه الشرب الشركة القابضة كله

قام بالمشاركة وابداء الرأى لهذا البرنامج:

قام بالتنسيق الفنى والإخراج لهذا الإصدار:

كيميائي/ محمود جمعه

الإدارة العامة للمسار الوظيفي- الشركة القابضة

للاقتراحات والشكاوى قم بمسح الصورة (QR)



